

Microbiological Culture Solutions



Nigeria Headquarter Global warehouse Training center Global subsidiary Manufacturing base

Haier Biomedical UK Ltd.















CONTENT

| • | CO ₂ Incubator | 03 |
|---|---|----|
| • | Standard Incubator | 15 |
| • | Climate Chamber | 19 |
| | Cell Culture Scene Consumables Management | 25 |
| • | Pipette Series | 30 |
| • | Centrifuge Tube Series | 32 |


CO₂ Incubator

Haier Biomedical IoT enabled CO₂ incubator with 180°C dry heat sterilisation provides a safe and secure reproducible growth environment for cell cultures.



IR Sensitive Control of CO₂ Concentration

The new IR sensor with high temperature resistance of 190°C is based on the NDIR measurement principle and uses a silicon MEMS transmitter to replace the traditional light source. It can withstand more than 300 dry heat sterilization cycles with a service life of up to 15 years and control accuracy of $\pm 0.1\%$. German IR infrared sensing technology, zero drift, without need for calibration, drift less than 0.3% within 2 years





7-inch Touchscreen

Displays CO2 concentration and temperature data in real time. 15 years of data can be exported via USB



Inner Door

The door ensures the inside of the cabinet is sealed

Outer Door

The heated outer door prevents the condensation of the inner door

Internal Partition

Safety anti-slip design of pull out shelves



Adjustable Feet It can be double stacked


304 Stainless Interior

180°C Dry-heat Sterilization

All internal components do not need to be disassembled and do not need separate autoclave sterilization to prevent secondary pollution. Cleaning consumables are not needed, one-button sterilization. German INFRARED CO2 sensor, NDIR light source technology drift < 0.3% within two years. The unit can withstand sterilization at 180oc with no disassembly and no manual calibration

Precise and Accurate Temperature Control

Controls the temperature precisely, within ±0.1°C, with six-sided heating based on the fuzzy PID control principle, to provide a stable temperature to ensure the normal growth of cells throughout their life cycle.



Temperature°C 37.1 37.0 36.9


Uniformity of 27 measuring points $<\pm 0.3^{\circ}C$

Central consistency point <±0.1°C

Precise CO₂ Concentration Using New IR Sensor Control Technology

Haier Biomedical's new IR Sensor technology uses NDIR measurement principles and withstands high temperatures of 190°C. The silicon MEMS transmitter can carry out more than 300 dry heat sterilization cycles to extend the service life to 15 years. Built-in temperature and humidity compensation technology reduces the impact of changes in humidity and temperature without the need for calibration after the high temperature sterilization. Five point calibration yields a higher measuring accuracy, sensitivity with less drift.







Silicon-based mems transmitter

Fast Environment Recovery for Optimal Cell Growth

Adopting active air flow control technology, and based on the fuzzy PID control principle, the parameters can be restored without overshoot. After opening the door for 30 seconds, the temperature and CO₂ concentration can be quickly restored within 4 minutes. Even if multiple users share a CO2 incubator and frequently open and close the door, the stability and uniformity of the incubator can be ensured.








Illustration of purified airflow

Temperature recovery curve (door open for 30s)

CO₂ concentration recovery curve (door open for 30s)

180°C Dry-Heat Sterilization Technology Minimises Contamination

Easy and effective sterilization of microorganisms including bacteria, fungi and microplasma with strong resistance, at 180°C high temperatures without the need for consumables. Simply press the "sterilization key" to activate and complete the sterilization process automatically in just 12 hours.

Delivers sterility level within the chamber of all surfaces to meet WS/T367-2012 standards.

All components are sterilized during the process, there is no need to dissemble internal components (including CO₂ sensors) and decontaminate separately, thus avoiding secondary pollution.

Temperature°C 200 Warm up 180 160 Sterilization 140 120 Cool down 100 80 60 40 20 10 11 12

High Efficiency Microbial Filter



The CO2 inlet is equipped with a high-efficiency microbial filter, with 99.99% filtration efficiency for particles larger than or equal to $0.2\mu m$ in diameter. It can effectively filter bacteria and dust particles in the CO2 gas line to ensure the safety of experimental results.

Easy to Clean Interior





The working chamber is plasma electro polished, stamped stainless steel with wide-arc, laser welded corners. Bracketless shelving design ensures that it is quick and easy to clean.

Interactive Intelligent Display with Easy Touch Operation

Touch-sensitive screen with rapid sensing even in rubber gloves. Green indicates normal operational parameters, while a red warning display indicates abnormal, making it easy to view data at a glance. A red warning display and audible buzzer will alarm when water level is low.



Home screen red warning.



Announcement function designed for multiple persons to use the same incubator making it clear to all users on important matters.



Real-time display of operation data & real-time display of temperature, for CO₂ concentration and O₂ concentration, and the data during the culture cycle can be viewed at any time.



Operation mode clear management authority: three-levels of authority to ensure the security of data.

Real-time Monitoring



An IoT module with multi-screen interface provides real-time uploadset parameters, operation parameters, operation curves, records and event records through the IoT cloud platform. The operation ofincubator can be monitored at anytime and anywhere throughmobile APP or computer terminal. Alarm function, and servicefunction are available through an one button touch.

Anti-Condensation Heating System to Reduce Pollution Risk

The door on the CO₂ incubator radiates heat to the inner glass door, effectively preventing the glass door from forming condensation. The possibility of microbial contamination caused by the condensate water is eliminated.

Intelligent Control of Circulating Air Maintains Uniformity

Automatically adjusts the circulation of the air flow, optimising the air flow to avoid air volatilization of samples and ensuring proper uniformity throughout the chamber.

Comprehensive Safety Alarm System

The system ensures the safety of experiments and processes by utilizing an independent temperature alarm system, including a sound light and remote reminder. Other alarms include CO₂ concentration, door ajar and water shortage.

CO₂ Incubator


Innovative and User-friendly Design with Attention to Detail



Safe anti-slip design with pull out shelves.



Drainage design




Active heat pipe condensation technology with any condensation directly returning to the reservoir.



Data traceable for 15 years with large storage capacity and data exportable through USB.

Accurate oxygen control (optional) for HCP-168



Three or six internal doors are available to reduce gas consumption

 O_2 concentration can be controlled within the range of 1-21% After opening the door for 30 seconds, the O_2 concentration can recover to 5% in only 8 minutes and 1% in 18 minutes High precision zirconia O_2 sensor, oxygen control accuracy 0.1% Advanced and reliable gas solenoid valve, low noise

The Quality of ISO Class 5 Clean Room Can Ensure a Better Cell Growth Environment



The optional HEPA high-efficiency filtration system combined with the unique air duct circulation design can continuously filter pollutants (biological pollutants and suspended particles) in the cabinet, ensuring that the incubator can reach the ISO class 5 clean room within 5 minutes after the external door is closed, which is equivalent to the class 100 environment of the 209 E standard of the united states

Optional Accessories

| Name | Material Description |
|-------------------------|--|
| Oxygen module | Zirconia O₂ sensor, control accuracy: 0.1%; control range: 1-21% |
| 3 Inner door | Reduce the temperature, humidity and carbon dioxide concentration in the box after opening the door, and minimize the mutual influence of multiple cultures |
| 6 Inner door | Reduce the temperature, humidity and carbon dioxide concentration in the box after opening the door, and minimize the mutual influence of multiple cultures |
| Water tray | Provides different bottom humidification methods |
| Roller base | Easy to move, prevent the ground bacteria contamination |
| HEPA filter | Ensure the cleanliness of the cabinet, suitable for users who open and close the door frequently; After opening the door for 30 seconds, the air inside the cabinet can be passed through HEPA filters within 5 minutes and reach ISO 5 clean room quality |
| Pressure reducing valve | Suitable for users with cylinder gas supply |
| Partition | Increase the number of samples cultured |

Specifications


| | Model | | HCP-80 HCP-80B | HCP-168 HCP-168B | HCP-258 HCP-258 | | |
|-----------------|---|-------|----------------------|-----------------------------------|---------------------|--|--|
| Гуре | | | Air Jacket | Air Jacket | Air Jacket | | |
| | Chamber Volume (L) | | 80 | 170 | 258 | | |
| | Interior Chamber | | 304 Stainless Steel | | | | |
| Construction | Exterior Chamber | | | Cold-rolled steel powder coated | d | | |
| | Access Port | | / | 35mm Diameter | 35mm Diameter | | |
| | Data Outputs | | Remote | Alarm Contacts, USB, and Option | nal 4-20mA | | |
| | | kg | 75/100 | 95/130 | 110/160 | | |
| | Net/Gross Weight (approx) | lbs | 165/220 | 209.4/286.6 | 243/353 | | |
| | | mm | 400*420*490 | 490*560*650 | 570*610*745 | | |
| _ | Interior Dimensions (W*D*H) | in | 15.7*16.5*19.3 | 19.3*22*25.6 | 22.4*24.0*29.3 | | |
| Dimensions | | mm | 625*684*735 | 714*812*887 | 794*867*985 | | |
| | Exterior Dimensions (W*D*H) | in | 24.6*26.9*28.5 | 28.1*32*34.9 | 31.3*34.1*38.8 | | |
| | | - | 765*710*930 | 890*800*1085 | 950*880*1145 | | |
| | Packing Dimensions (W*D*H) | mm | | | | | |
| | Discossions (M*D) | in | 30*28*36.6 | 35*31.5*43 | 37.4*35*45 | | |
| | Dimensions (W*D) | mm | 380*300 | 473*434 | 550*484 | | |
| Shelves | Number Standard/Maximum | Ι. | 3/8 | 3/11 | 3/13 | | |
| | Max.Load Per Shelf/Total Load | kg | 15/45 | 15/45 | 15/45 | | |
| | Construction | | | Perforated, Adjustable | | | |
| Electrical | Rated Voltage Power Supply (V/Hz) | | 220-240~50/60 115/60 | 220-240~50/60 115/60 | 220-240~50/60 115/6 | | |
| | Steri-Run Consumption (kw) | | 0.85/0.75 | 1.3/1.1 | 1.35/1.2 | | |
| Control | Controller | | Microprocessor | Microprocessor | Microprocessor | | |
| COILLOI | Display | | 7 "LCD Screen | 7 "LCD Screen | 7 "LCD Screen | | |
| | Control Accuracy | | 0.1% | 0.1% | 0.1% | | |
| | Range | | 0-20% | 0-20% | 0-20% | | |
| | Alarm Range | | ±0.5% | ±0.5% | ±0.5% | | |
| | Inlet Pressure | | | 12-17Psi (0.8-1.2 Bar) | | | |
| CO ₂ | Gas Purity | | ≥99.5% | | | | |
| | Sensor | | IR | IR | IR | | |
| | Recovery Time at 5vol%/CO ₂ | | 4 | 4 | 4 | | |
| | for a 30 Second Door Opening* (min) | | -0.2 | | | | |
| | CO ₂ Inlet Filter (µm) | | <0.2 Y | <0.2 Y | <0.2 Y | | |
| | High/Low Temperature | | | | | | |
| | Remote Alarm | | Y | Y | Y | | |
| Alarms | Excessive CO ₂ Concentration | | Y | Y | Y | | |
| | Water Shortage | | Y | Y | Y | | |
| | Sensor Error | | Y | Y | Y | | |
| | Door Ajar | | Y | Y | Y | | |
| | Control Accuracy (°C) | | 0.1 | 0.1 | 0.1 | | |
| | Range | | | Ambient temperature+3-55°C | | | |
| Temperature | Uniformity (°C) | | ±0.3 | ±0.3 | ±0.3 | | |
| Parameter | Ambient Range (°C) | | 18-32 | 18-32 | 18-32 | | |
| | Sensor Recovery Time at 37°C | | 2PT1000 | 2PT1000 | 2PT1000 | | |
| | for a 30 Second Door Opening* | (min) | 4 | 4 | 4 | | |
| Sterilization | Cycle Temperature | | 18 | 30°C on internal Surfaces and She | elves | | |
| Cycle | Cycle Duration | | Under 12 Hours | Under 12 Hours | Under 12 Hours | | |
| Llumidity | RH (Relative Humidity) | | Setting 37°C ≥90% | Setting 37°C ≥90% | Setting 37°C ≥90% | | |
| Humidity | Humidity Reservoir | | Max. 1.3L/Min 0.5L | Max. 3L/Min 0.5L | Max. 3.6L/Min 0.5L | | |
| | Hepa Filter | | Optional | Optional | Optional | | |
| | Pressure Reducing Valve | | Optional | Optional | Optional | | |
| Accessories | RS485 | | ү | Y | Y | | |
| | 4-20mA | | Optional | Optional | Optional | | |
| | Cylinder Switching Devicey | | Optional | Optional | Optional | | |
| Certification | Symbol Symbol in ig Devicey | | CE UL | CE UL | CE UL | | |

Product appearance and specifications are subject to change without notice




Precise and Accurate Temperature Control

Controls the temperature precisely, within ± 0.1 °C, with six-sided heating based on the fuzzy PID control principle, to provide a stable temperature to ensure the normal growth of cells throughout their life cycle.



Precise CO₂ Concentration Using New IR Sensor Control Technology


Haier Biomedical's new IR Sensor technology uses NDIR measurement principles and with stands high temperature of 100° C. The silicon MEMS transmitter can carry out more than 300 dry heat sterilization cycles to extend the service life to 15 years. Built-in temperature and humidity compensation technology reduce the impact of changes in humidity and temperature without the need for calibration after the high temperature sterilization. Five points calibration yields a higher measuring accuracy, sensitivity with less drift (less than 3% within 2 years).



^{*}The equipment is tested by Haier in a controlled environment. Haier does not guarantee that the results of field tests under different conditions will be consistent. The test model is HCP-168E

Fast Environment Recovery for Optimal Cell Growth

Adopting active air flow control technology, based on the fuzzy PID control principle, the parameters can be restored without overshoot. After opening the door for 30 seconds, the temperature and CO₂ concentration can be quickly restored within 4 minutes. Even if multiple users share a CO₂ incubator and frequently open and close the door, the stability and uniformity of the incubator can be ensured.




90°C Moist Heat Sterilization Technology

Effective sterilization of microorganisms including bacillus and spores with strong resistance, at 90° C under moist heat, without the need for consumables. Simply press the "sterilization button", to activate and complete the sterilization process automatically in 14 hours.

Delivers sterility level within the chamber of all surfaces to meet WS/T367-2012 standards.

All components are sterilized during the process, there is no need to dissemble internal components (including CO₂ sensors) and decontaminate separately, thus avoiding secondary pollution.



Sterilization Temperature Profile

Forty-seven points were tested in the working chamber, including glass inner doors and partitions. All regions reached 90°C and maintained for 9 hours.

^{*}The equipment is tested by Haier in a controlled environment. Haier does not guarantee that the results of field tests under different conditions will be consistent.

The test model is HCP-168E

CO₂ Incubator

Door Switch

When the door opens, heating, air intake and fan automatically stop to minimize the risk of cross contamination

Co₂ Sensor

- \bullet The new IR sensor with high temperature resistance of 100 °C , can withstand more than 300 high heat sterilization cycles
- Based on the NDIR measurement principle and uses a silicon MEMS transmitter to replace the traditional light source
- Zero drift and without need for calibration



Partition

- Anti-slip design
- High levelness ensures uniform growth of adherent cells
- Mirror stainless steel to ensure high surface cleanliness, easy to clean



Air Flow System

The air flow circulation ensures proper uniformity throughout the chamber

Integrated Liner

Integral design, large arc design, easy to clean $\,$



Air Jacketed With Six-sides Heating Design

- \bullet Fast temperature recovery and superior temperature uniformity
- High temperature sterilization can ensure that the temperature of each surface can reach 90°C

Inner Door

- •Tempered glass provides easy observation of sample growth
- •Three/six inner doors optional

Operation Panel

- 4-inch LCD screen, vivid display and easy operation
- Abnormal operation sound and light alarm to ensure sample safety
- Running data can be traced, large capacity storage, data can be exported through USB



Test Hole

Providing access for convenient measurement of internal statistics



Outer Door

- Prevents the condensation of the inner door
- Left/right hand door optional

Inner and Outer Door Seal

- Silicone material, prevent aging after heating
- Close the inner cavity to ensure the cleanliness and uniformity of the inner chamber


Bottom Reservoir Humidification

- Reservoir humidification method, no water tray, easy to clean, avoid breeding bacteria
- Large evaporation area and fast humidity recovery



Easy to Clean Interior

The working chamber is plasma electro polished, stamped stainless steel with wide-arc, laser welded corners. Bracketless shelving design ensures that it is quick and easy to clean.





Innovative Design with Attention to Detail






Safe anti-slip design of pull-out shelves.

Data traceable for 15 years with large storage capacity and data exportable through USB.

Reservoir Humidification Without Condensation

Active heat pipe condensation technology with condensate water directly returns to the reservoir, to ensure no condensation.



Optional Accessories

| Name | Material Description |
|-------------------------|---|
| Oxygen module | Zirconia O₂ sensor, control accuracy: 0.1%; control range: 1-21% |
| 3 Inner door | Reduce the temperature, humidity and carbon dioxide concentration in the box after opening the door, and minimize the mutual influence of multiple cultures |
| 6 Inner door | Reduce the temperature, humidity and carbon dioxide concentration in the box after opening the door, and minimize the mutual influence of multiple cultures |
| Water tray | Provides different bottom humidification methods |
| Roller base | Easy to move, prevent the ground bacteria contamination |
| Stacking frame | Stacking the two incubators makes the fixation firmer |
| HEPA filter | Ensure the cleanliness of the cabinet, suitable for users who open and close the door frequently; After opening the door for 30 seconds, the air inside the cabinet can be passed through HEPA filters within 5 minutes andreach ISO 5 clean room quality |
| Pressure reducing valve | Suitable for users with cylinder gas supply |
| Partition | Increase the number of samples cultured |


Specifications

| | Model | Power Supply (V/Hz) | Volume (L) | Exterior Dimensions (W*D*H)(mm) | Interior Dimensions (W*D*H)(mm) | Shelf Dimensions (W*D*H)(mm) | Standard Configuration of Shelves No./Maximum |
|---|----------|------------------------|------------|------------------------------------|------------------------------------|---------------------------------|--|
| | HCP-80E | 220-240~50/60 | 80 | 625*684*735 | 400*420*490 | 380*300 | 3/8 |
| ŀ | HCP-168E | 220-240~50/60 | 170 | 714*812*887 | 490*560*650 | 473/434 | 3/11 |
| H | HCP-258E | 220-240~50/60 | 258 | 794*867*985 | 570*610*745 | 550*484 | 3/13 |

Standard Incubator

Scope of Application

The solution is widely used in bacteria, fungi and other microorganisms culture; as well as enzyme digestion reaction, ligation reaction, embedded incubation and other related constant temperature experiments.







Embedded incubation

Bacter

Fungus



Product Advantages



Personalized interface, easy to link

Equipped with USB and RS485 interfaces to meet the different needs of users to transfer data



Multiple protection benefits for increased security

Overheat protection (OPT), over current protection (FU), sensor error detection, independent temperature limit, compliance with DIN 12880 requirements and EU 3.1 safety level. Sound, light and remote alarms (optional) which guarantee experiment safety. Multiple alarms, such as over temperature alarm, high and low temperature alarm, door ajar, and sensor error alarm.



Data traceability

Data traceable up to 15 years with base storage 8GB (64GB optional) and data exportable through USB



High thermal insulation performance, energy saving and environmental protection

The unit is manufactured with aluminum foil insulation cotton, which improves the overall insulation performance and reduces energy consumption, lowering costs while also being environmentally friendly.

Fuzzy PID Control Technology







ASTM standard, 12 points testing


Based on PID control principle, manufactured with U-shaped 3-sided heating to achieve superior temperature control and uniformity control.


Rapid Recovery After Door Open

Natural Convection Forced Convection

Rapid warming: the temperature inside the unit quickly recovers after opening the door to reduce the influence of temperature fluctuation on the sample.

The temperature rise curve to 37°C after opening the door for 30 sec at 22°C ambient temperature





Convenient and Intelligent Management at a Glance



7-inch touchscreen, easy to operate and sensitive, it can respond quickly even when wearing rubber gloves.



Real-time display of temperature data, one-touch to review previous data.



Records abnormal information in real time, eliminating any hidden abnormalities which ensures the culturing is more secure.



Multiple operating modes.




The program can be edited and set at any number of segments to meet the needs of various detection tests.


Standard Incubator

Precise Temperature Control, Energy-efficient and Environment-friendly

An energy-efficient model with superior control and heating mechanisms, high-quality insulation material and cabinet structure to ensure heating requirements are met while keeping power consumption to a minimum.



HZP-168



Optional IoT Technology for Real-time Remote Monitoring



Through the mobile app, the status of the incubator can be checked in real time, and information such as temperature alarm, sensor error alarm and door ajar can be controlled with one button, which provides more security for the experiment process.

Pictures in Details



Seamless, curved internal chamber for easy cleaning and decontamination.



Standard independent intelligent temperature safety controller to ensure experimental safety; RS485 achieves seamless IoT data connection.

Product Parameters

| Model | Product Series | Capacity (L) | Power (W) | Exterior Dimensions (W*D*H)(mm) | Interior Dimensions (W*D*H)(mm) | Packing Dimensions (W*D*H)(mm) |
|---------|--------------------|-----------------|--------------|------------------------------------|------------------------------------|-----------------------------------|
| HZP-168 | Natural convection | 168 | 640 | 650*782*1028 | 490*550*626 | 795*835*1220 |
| HFP-80 | Forced convection | 80 | 510 | 560*662*870 | 400*400*480 | 710*760*1070 |

| Shelves (Standard) | Temperature Control Range | Temperature Uniformity | Temperature Fluctuation (°C) | Temperature Control Precision (°C) | Recovery Time after 30 sec Door Opening (min) |
|-----------------------|---------------------------|------------------------|------------------------------|------------------------------------|--|
| 2/17 | RT+5~105°C | ±0.5°C at 37°C | ±0.1 | ±0.1 | 5 |
| 2/12 | RT+5~105°C | ±0.3°C at 37°C | ±0.1 | ±0.1 | 2.5 |

Product appearance and specifications are subject to change without notice

Animal and plant tissue culture, drug stability tests, cosmetic stability tests, food shelf life tests, electronic components aging tests, packaging material stability tests.



HHS-256/756/506

Product Advantages



Silent

Semiconductor technology ensures low vibration and noise output with no pollution into the environment



Precise control

Accurate temperature and humidity control, long-term stability, 40 °C temperature uniformity ±0.5°C and central temperature fluctuation ±0.2°C, 75% humidity fluctuation ±1%



Water-saving

Intelligent control of PTC humidification, daily water consumption of 120-320ml, no need to recycle waste water, saving space



Power saving

Semiconductor technology means the daily power consumption is as low as 5kWh; 90% more energy efficient than compressor technologies

Product Features



Multiple protection protocols - equipped with delay start, high/low temperature and light intensity protection in line with DIN12880 requirement for over/under temperature protection



Expandable large capacity data storage, the touchscreen memory can be expanded to 64GB, storing up to 15 years data which can be exported via a USB



High precision temperature sensor, dual PT1000 sensors for more accurate temperature control



An access port with a diameter of 35mm on the left side of the cabinet to facilitate independent testing of temperature and humidity



Optional electromagnetic lock, suitable for multiple users with independent management for safety



High insulating performance polyethene foam provides excellent insulation and stable cabinet temperatures reducing energy consumption



Polyurethane foam insulation provides excellent thermal insulation, reducing energy consumption



High precision capacitive humidity sensor




Microprocessor control system


- PID control principle, 10-inch touch screen, temperature control precision 0.1°C, humidity control precision 0.1%, temperature range 5-70°C, humidity range 10%-90%
- USB, RS485 interface as standard
- Temperature alarm, humidity alarm, door alarm, sensor alarm and water shortage warning
- Display temperature, humidity and ambient temperature; users can query the historical curve

International Quality Assurance

Accurate Temperature Control



DIN12880 standard 27 test points



27 test points temperature uniformity is ± 0.5 °C (ambient temperature: 22°C, setting temperature: 40°C)


Haier Biomedical

Product Parts



Climate Chamber


Accurate Humidity Control



A: Stable temperature and relative humidity control range


B: Ambient temperature: lower than 22 ° C, ambient humidity: lower than 40%

■: ICH standard test points





 $Humidity fluctuation \pm 1\% \\ (setting temperature: 40°C, setting humidity: 75\%. \\ ambient temperature: 22°C, ambient humidity: 40%)$

Internal Cabinet Environment Quick Recovery System



Professional air duct design, ensuring temperature and humidity uniformity





The time required for humidity recovery to 75% after opening door for 30s is <14 minutes

Intelligent Management

• Convenient and intelligent management to improve working efficiency



The intelligent 10-inch touchscreen controller is easy to operate and sensitive to touch, even in rubber gloves. The PID control algorithm ensures the accuracy of temperature control



Data and multi-user authority management and permissions conforms to FDA 21 CFR Part 11



Unlimited programs with infinite humidity and temperature settings to allow users to customise to their needs

• High quality manufacture and reliable operation



Capacitive humidity sensor, long-term operating reliability

- Interference-free humidity data collection.
- Long-term reliability without the need for calibration.
- High precision $\pm 0.1\%$.
- Anti-condensation design for more accurate humidity monitoring



High precision temperature sensor, accurate and reliable

- Adopts PT1000 temperature sensors for accurate, stable and repeatable measurement without deviation.
- Dual sensors further improve accuracy.



Semi-conductor cooling, superior energy-saving and mute effect

 Semiconductor thermocouple consists of N-shape semiconductor and P-shape semiconductor.

• Intelligent control, ensures temperature and humidity accuracy



Intelligent control PTC humidification, energy-saving and water-saving $\,$

The temperature and purity of vapour are accurately controlled by the intelligent water supply system and ceramic high-temperature heating apparatus.



Intelligent dehumidification, accurate humidity control

Semi-conductor intelligent dehumidification system accurately controls heating and cooling, matching with humidity control.

Specifications

| | Model | | HHS-256 | HHS-506 | HHS-756 |
|--------------|--|-------|---------------------------------|---------------------------------|---------------------------------|
| | Chamber Volume (L) | | 256L | 506L | 756L |
| 0 | Interior Chamber | | stainless steel | stainless steel | stainless steel |
| Construction | Exterior Chamber | | Galvanized Sheet Powder Coating | Galvanized sheet powder coating | Galvanized sheet powder coating |
| | Access Port | | 35mm Diameter | 35mm Diameter | 35mm Diameter |
| | Net/Gross Weight | kg | 175/188 | 225/260 | 280/328 |
| Dimensions | Interior Dimensions (W*D*H) mm | | 650*570*700 | 740*570*1200 | 1100*570*1200 |
| DIFFICUSIONS | Exterior Dimensions (W*D*H) | mm | 835*905*1190 | 930*905*1690 | 1290*905*1690 |
| | Packing Dimensions (W*D*H) | mm | 1030*955*1280 | 1110*955*1780 | 1380*955*1780 |
| | Dimension/mm(W*D) | | 597*531 | 687*531 | 1048*531 |
| Shelves | Standard Qty / Max Qty | | 2/16 | 2/31 | 2/31 |
| SHEIVES | Max Weight Per Shelf | kg | 20 | 20 | 20 |
| | Structure | | Slide rail, adjustable | Slide rail, adjustable | Slide rail, adjustable |
| | Voltage / Frequency (V/Hz) | | 220-240~50/60 | 220-240~50/60 | 220-240~50/60 |
| Electrical | Power (W) | | 750 | 1100 | 1760 |
| | Day Consumption at 25°C & 40% RH (kw·h) | | 4.6 | 5.4 | 5.6 |
| Control | Controller | | The microprocessor | The microprocessor | The microprocessor |
| | Display | | 10 "smart LCD screen | 10 "smart LCD screen | 10 "smart LCD screen |
| | The Set Range (°C) | | 5~70 | 5~70 | 5~70 |
| | Control Precision (°C) | | ±0.1 | ±0.1 | ±0.1 |
| The | Temperature Uniformity at 25 °C | | ±0.2 | ±0.2 | ±0.2 |
| Temperature | Temperature Fluctuation at 25 °C | | ±0.1 | ±0.1 | ±0.1 |
| Parameter | The Sensor | | PT1000 | PT1000 | PT1000 |
| | Rate of Temperature Rise (°C / min) | | 1 | 0.8 | 0.6 |
| | 30 Seconds Recovery Time After Door Opening at 40°C (min) | | 3 | 3.8 | 5 |
| | Humidity Setting Range (% RH) | | 10~90 | 10~90 | 10~90 |
| Humidity | Humidity Setting Accuracy (% RH) | | 0.1 | 0.1 | 0.1 |
| Parameter | Humidity Fluctuation at 25 °C & 40% RH | (%RH) | ±0.5 | ±0.5 | ±0.5 |
| | Daily Water Consumption (ml) | | 120 | 240 | 320 |
| Ontional | Electromagnetic lock (password) | | Υ | Y | Υ |
| Optional | Printer | | Y | Y | Y |
| | Remote Alarm Interface | | Y | Y | Υ |
| Standard | RS485 | | Y | Y | Y |
| | Water Level Alarm | | Y | Y | Y |

Product appearance and specifications are subject to change without notice $\,$



Scope of Application

In order to fully satisfy the growth, survival and reproduction of cells in vitro without pollution, the selection of consumables is very important. Haier Biomedical culture consumables meet the high-quality, full range of cell culture scenarios, is the ideal choice for cell culture, and comprehensively contributes to the development of life science. The raw materials of the consumables products meet the USP VI standard to ensure the qualification of the production of raw materials. The equipment from Germany and Japan is used to achieve a variety of non-invasive production processes of consumables. The whole process of 100,000 GMP purification workshop production, cobalt-60 irradiation sterilization, fully ensure the sterility and uniformity of products, all-round guarantee the cell growth process is safe and clean.

Cell Culture Flask



- High quality polystyrene material, 100,000 clean workshop manufacturing
- TC treated and non-TC treated specifications can be used for cell adhesion culture and cell suspension culture
- Large mouth design for easy removal of cell suspension or digestive cells
- Two closure type to choose: One type is with septum cap, which has a 0.22µm membrane that allows sterile air exchange or a leak resistant seal, the other is with screw cap
- Product irradiation sterilization. RNase, DNase, Pyrogen and Endotoxin Free
- Scales and marking areas on both sides of the bottle

| CAT.NO. | Description | Unit Quantity | Packing Specification |
|---------|---|---------------|------------------------|
| 633017 | 25cm², sealed cap, wall culture, TC treatment | 200 | 10pcs/bag, 20bags/case |
| 634018 | 25cm², vented cap, wall culture, TC treatment | 200 | 10pcs/bag, 20bags/case |
| 633019 | 75cm², sealed cap, wall culture, TC treatment | 100 | 5pcs/bag, 20bags/case |
| 634020 | 75cm², vented cap, wall culture, TC treatment | 100 | 5pcs/bag, 20bags/case |
| 633021 | 175cm², sealed cap, wall culture, TC treatment | 40 | 5pcsbag, 8bags/case |
| 634022 | 175cm², vented cap, wall culture, TC treatment | 40 | 5pcsbag, 8bags/case |
| 633023 | 225cm², sealed cap, wall culture, TC treatment | 25 | 5pcs/bag, 5bags/case |
| 634024 | 225cm², vented cap, wall culture, TC treatment | 25 | 5pcs/bag, 5bags/case |
| 631025 | 25cm², sealed cap, suspension culture, untreated | 200 | 10pcs/bag, 20bags/case |
| 632026 | 25cm², vented cap, suspension culture, untreated | 200 | 10pcs/bag, 20bags/case |
| 631027 | 75cm², sealed cap, suspension culture, untreated | 100 | 5pcs/bag, 20bags/case |
| 632028 | 75cm², vented cap, suspension culture, untreated | 100 | 5pcs/bag, 20bags/case |
| 631029 | 175cm², sealed cap, suspension culture, untreated | 40 | 5pcs/bag, 8bags/case |
| 632030 | 175cm², vented cap, suspension culture, untreated | 40 | 5pcs/bag, 8bags/case |
| 631031 | 225cm², sealed cap, suspension culture, untreated | 25 | 5pcs/bag, 5bags/case |
| 632032 | 225cm², vented cap, suspension culture, untreated | 25 | 5pcs/bag, 5bags/case |

Erlenmeyer Flask



Product Features

- Two type of materials are available: PC/PETG with baffled bottom/flat bottom
- Two closure type to choose: One type is with septum cap, which has a 0.22µm membrane that allows sterile air exchange or a leak resistant seal, the other is with screw cap
- Single independent vacuum packaging, easy to use, avoid cross contamination
- Product irradiation sterilization. RNase, DNase, Pyrogen and Endotoxin Free
- The scale on the outer wall of the bottle is clear and precise, and the inner wall is smooth
- Both suspension cells and adherent cells can be cultured

| CAT.NO. | Description | Unit Quantity | Packing Specification |
|---------|--------------------------|---------------|-----------------------|
| 641049 | 125ml, PETG, sealed cap | 24 | 1pcs/bag, 24bags/case |
| 641050 | 250ml, PETG, sealed cap | 12 | 1pcs/bag, 12bags/case |
| 641051 | 500ml, PETG, sealed cap | 12 | 1pcs/bag, 12bags/case |
| 641052 | 1000ml, PETG, sealed cap | 12 | 1pcs/bag, 12bags/case |
| 642053 | 125ml, PETG, vented cap | 24 | 1pcs/bag, 24bags/case |
| 642054 | 250ml, PETG, vented cap | 12 | 1pcs/bag, 12bags/case |
| 642055 | 500ml, PETG, vented cap | 12 | 1pcs/bag, 12bags/case |
| 642056 | 1000ml, PETG, vented cap | 12 | 1pcs/bag, 12bags/case |
| 643057 | 125ml, PC, sealed cap | 24 | 1pcs/bag, 24bags/case |
| 643058 | 250ml, PC, sealed cap | 12 | 1pcs/bag, 12bags/case |
| 643059 | 500ml, PC, sealed cap | 12 | 1pcs/bag, 12bags/case |
| 643060 | 1000ml, PC, sealed cap | 12 | 1pcs/bag, 12bags/case |
| 644061 | 125ml, PC, vented cap | 24 | 1pcs/bag, 24bags/case |
| 644062 | 250ml, PC, vented cap | 12 | 1pcs/bag, 12bags/case |
| 644063 | 500ml, PC, vented cap | 12 | 1pcs/bag, 12bags/case |
| 644064 | 1000ml, PC, vented cap | 12 | 1pcs/bag, 12bags/case |
| 645065 | 125ml, PC, sealed cap | 24 | 1pcs/bag, 24bags/case |
| 645066 | 250ml, PC, sealed cap | 12 | 1pcs/bag, 12bags/case |
| 645067 | 500ml, PC, sealed cap | 12 | 1pcs/bag, 12bags/case |
| 645068 | 1000ml, PC, sealed cap | 12 | 1pcs/bag, 12bags/case |
| 646069 | 125ml, PC, vented cap | 24 | 1pcs/bag, 24bags/case |
| 646070 | 250ml, PC, vented cap | 12 | 1pcs/bag, 12bags/case |
| 646071 | 500ml, PC, vented cap | 12 | 1pcs/bag, 12bags/case |
| 646072 | 1000ml, PC, vented cap | 12 | 1pcs/bag, 12bags/case |

Cell Culture Plate



- A wide range of sizes, materials and shapes to suit different experimental needs of the laboratory
- There are TC treated and non-TC treated options which is suitable for cell suspension culture and cell adhesion culture
- Sterile, no enzyme, no heat source and endotoxin, individually packed to ensure the stability of the experimental results
- High transparency polystyrene material can be used for microscope observation
- Raised edge on wells to prevent cross contamination
- The left side and the upper side have alphanumeric labels, easy to distinguish and identify, easy to record
- Suitable for most plate readers and automated equipment

| CAT.NO. | Description | Unit Quantity | Packing Specification |
|---------|--|---------------|------------------------|
| 612001 | 96-well, blister box/F-bottom (flat bottom), TC treatment, wall culture | 100 | 1pcs/bag, 100bags/case |
| 612002 | 48-well, blister box/F-bottom (flat bottom), TC treatment, wall culture | 100 | 1pcs/bag, 100bags/case |
| 612003 | 24-well, blister box/F-bottom (flat bottom), TC treatment, wall culture | 100 | 1pcs/bag, 100bags/case |
| 612004 | 12-well, blister box/F-bottom (flat bottom), TC treatment, wall culture | 100 | 1pcs/bag, 100bags/case |
| 612005 | 6-well, blister box/F-bottom (flat bottom), TC treatment, wall culture | 100 | 1pcs/bag, 100bags/case |
| 611006 | 96-well, blister box/F-bottom (flat bottom), untreated, suspension culture | 100 | 1pcs/bag, 100bags/case |
| 611007 | 48-well, blister box/F-bottom (flat bottom), untreated, suspension culture | 100 | 1pcs/bag, 100bags/case |
| 611008 | 24-well, blister box/F-bottom (flat bottom), untreated, suspension culture | 100 | 1pcs/bag, 100bags/case |
| 611009 | 12-well, blister box/F-bottom (flat bottom), untreated, suspension culture | 100 | 1pcs/bag, 100bags/case |
| 611010 | 6-well, blister box/F-bottom (flat bottom), untreated, suspension culture | 100 | 1pcs/bag, 100bags/case |

Cell Culture Dish







Product Features

- Assorted sizes, styles and surface finishes for different culture requirements
- Designed with raised edges, easy to stack and transport
- $\, \bullet \,$ Polystyrene material, surface treatment, suitable for cell adhesion growth
- Vents designed to maintain gas exchange under sterile conditions
- Sterile, Enzyme, Pyrogen and Endotoxin Free, individual packaging, to ensure the stability of the experimental results

| CAT.NO. | Description | Unit Quantity | Packing Specification |
|---------|---|---------------|------------------------|
| 622011 | 6cm, plastic sealed bag/easy grip, TC treated, wall culture | 500 | 10pcs/bag, 50bags/case |
| 622012 | 10cm, plastic sealed bag/easy grip, TC treated, wall culture | 300 | 10pcs/bag, 30bags/case |
| 622013 | 15cm, plastic sealed bag/easy grip, TC treatment, wall culture | 100 | 5pcs/bag, 20bags/case |
| 621014 | 6cm, plastic sealed bag/easy grip, untreated, suspension culture | 500 | 10pcs/bag, 50bags/case |
| 621015 | 10cm, plastic sealed bag/easy grip, untreated, suspension culture | 300 | 10pcs/bag, 30bags/case |
| 621016 | 15cm, plastic sealed bag/easy grip, untreated, suspension culture | 100 | 5pcs/bag, 20bags/case |


Serological Pipette



- Multiple volume sizes available, different color rings on the tube head for easy identification of different models
- Bi-directional scale design for easy reading of pipette volume, negative scale increases pipette capacity for larger volumes
- Clear and accurate scale with a precision of ±2% of total volume
- Filter cartridge prevents samples, aerosols and water vapor from entering the pipette to avoid cross-contamination
- Adaptable to a wide range of pipettes on the market with rubber adapters
- Multiple packaging options, easy-tear paper-plastic single packs as well as large packs for easy batch use and reduced packaging waste
- Irradiated, meet SAL 10⁻⁶, DNase/RNase free, no pyrogenic endotoxin

| CAT.NO. | Description | Unit Quantity | Packing Specification |
|---------|--|---------------|------------------------|
| 692109 | 1ml (yellow) disposable serum pipette, sterile | 1000 | 50pcs/bag, 20bags/box |
| 692110 | 2ml (green) disposable serum pipettes, sterile | 1000 | 50pcs/bag, 20bags/box |
| 692111 | 5ml (blue) disposable serum pipettes, sterile | 500 | 50pcs/bag, 10bags/case |
| 692112 | 10ml (orange) disposable serum pipettes, sterile | 500 | 50pcs/bag, 10bags/case |
| 692113 | 25ml (red) disposable serum pipettes, sterile | 250 | 25pcs/bag, 10bags/case |
| 692114 | 50ml (violet) disposable serum pipettes, sterile | 200 | 25pcs/bag, 8bags/case |
| 692115 | 100ml (black) disposable serum pipettes, sterile | 120 | 20pcs/bag, 6bags/case |

Pasteur Pipette



Product Features

- $\, \bullet \,$ Good fluidity of the tube wall, easy to control the pipetting volume during operation
- High transparency, clear scale, easy to observe the pipetting process in real time
- Made of high quality material, it can be used for pipetting liquids in all kinds of containers
- Individual packages, easy to use, well sealed, reproducible pipetting
- Each package is individually labeled with a lot number for quality tracking and traceability
- Enzymatic, Pyrogenic, Endotoxic, Cytotoxic and Hemolytic Free
- EO sterilized to ensure contamination-free experiments

| CAT.NO. | Description | Unit Quantity | Packing Specification |
|---------|--|---------------|--|
| 603001 | Total volume 1.5m. Each 0.1ml scale to 0.3ml, length 116mm, about 23 drops/ml, E.O. sterile, LDPE material | 1000 | Plastic pack, 100pcs/ box, 1,000pcs/case |
| 603002 | Total volume 3ml. Every 1/4ml scale to 1ml, length 140mm, E.O. sterile, LDPE material | 1000 | Plastic pack, 100pcs/ box, 1,000pcs/case |
| 603003 | Total volume 5ml. Every 1/4 ml scale to 1 ml, length 145 mm, E.O. sterile, LDPE material | 1000 | Plastic pack, 100pcs/box, 1,000pcs/case |
| 603004 | Total volume 5ml. Every 1/4ml scale to 1ml, length 150mm, E.O. sterile, LDPE material | 1000 | Plastic pack, 100pcs/ box, 1,000pcs/case |
| 603005 | Total volume 5ml. Every 1/2ml scale to 2ml, length 155mm, E.O. sterile, LDPE material | 1000 | Plastic pack, 100pcs/box, 1,000pcs/case |
| 603006 | Total volume 7ml. Every 1/2ml scale to 3ml, length 155mm, E.O. sterile, LDPE material | 1000 | Plastic pack, 100pcs/ box, 1,000pcs/case |
| 603007 | Total volume 7.5ml. Every 1/2ml scale to 3ml, length 148mm, approx. 23 drops/ml, E.O. sterile, LDPE | 1000 | Plastic pack, 100pcs/box, 1,000pcs/case |

Microcentrifuge Tubes



- Made of USP Class VI, medical-grade polypropylene (PP), free from release agent, plasticizer, bacteriostatic agent, and heavy metals
- High-clarity tube for convenient sample observation
- The unique cap design with flat surface for writing, delivers enhanced sealing and enables single-handed operation
- Temperature range: -80°C~120°C



Ultra-high clarity makes samples clearly visible
Universal design suitable for all major brands



Frosted writing area



Snap-on safety cap to ensure safe centrifugation to maximize operational safety during centrifugation

| CAT.NO. | Description | Unit Quantity | Packing Specification |
|---------|---|---------------|------------------------|
| 311001 | 0.6ml microcentrifuge tube, non-sterile | 10000 | 1000pcs/box, 10box/ctn |
| 311002 | 1.5ml microcentrifuge tube, non-sterile | 5000 | 500/box, 10boxes/case |
| 311004 | 2.0ml microcentrifuge tube, non-sterile | 5000 | 500/box, 10boxes/case |