

Laboratory Protection Clean Air Solutions

Your Haier Biomedical Partner



Qingdao Haier Biomedical Co., Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266111, P.R. China
E-mail: inquiry@haierbiomedical.com
Website: www.haiermedical.com
Dec, 2025



Haier Biomedical
International



Haier Biomedical
International



@haiermedicalint



Haier Biomedical
International



Haier Biomedical
International

Note: If a slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation of this brochure, please contact us for any further information if required.

CONTENTS

- » **Certifications, Quality Control, Patented Technologies** ----- 02
- » **Safe Admixture Solution for Medicine** ----- 03
- » **Differences between Biological Safety Cabinet and Laminar Flow Cabinet** --- 04
- » **Model Selection Guide for Biological Safety Cabinet** ----- 04
- » **Haier Biomedical Safety Cabinet, Type A2** ----- 05
 - NSF Series Biological Safety Cabinet ----- 05
 - N Series Class II Type A2 Biological Safety Cabinet ----- 08
 - Biological Safety Cabinet (Single/Double exhaust filter) ----- 13
 - Biological Safety Cabinet ----- 16
 - Mini Biological Safety Cabinet ----- 18
 - Classic Series Biological Safety Cabinet (Type A2) ----- 21
- » **Haier Biomedical Safety Cabinet, Type B2** ----- 24
- » **Animal Containment Workstation** ----- 26
- » **Haier Biomedical Laminar Flow Cabinet** ----- 27


Certifications, Quality Control, Patented Technologies

» Certifications


Safety certification to EN 61010
 EMC certification EN 61326
 Certified EN 12469 for Biological Safety Cabinets
 Certified CFDA YY-0569
 ISO 13485:2016 and ISO 9001:2015 Certified Company



Strict QC Tests and Pre-delivery Inspections



Down-flow velocity measuring



Inflow velocity measurement



UV lamp checks



Interior light checks



Workbench vibration checks



Filter integrity testing



Noise level checks

» Patented Technologies



LNS energy-saving mode (the fan will stop automatically once people leave for 15 minutes)



Intelligent constant air velocity



Pressure sensors monitoring service life of filters



UV lamp one-touch protocol



Prevent airflow from overflowing


Safe Admixture Solution for Medicine

>> Typical Application for PIVAS (Pharmacy Intravenous Admixture Service)


Haier Biomedical clean bench ensures a superior cleanliness environment while the technical specialists /medical staff perform the admixture of intravenous fluid for PIVAS.



>> Introduction to Safe System Solution for PIVAS




Complete Safe Solution for PIVAS



Differences between Biological Safety Cabinet and Laminar Flow Cabinet

>> Working Principles



>> Functions

| Product Category | Airflow Circulation | Applications | Air Supply Mode | Filter | Operator Protection | Sample Protection | Environment Protection |
|--|----------------------------------|---|--|-----------------|---------------------|-------------------|------------------------|
| Biological safety cabinet, Class II, Type A2 | 70% circulated, 30% discharged | Operation of pathogenic bacteria, mold, yeast and other hazardous samples | Negative pressure (Air pulled into cabinet) | High efficiency | ✓ | ✓ | ✓ |
| Biological safety cabinet, Class II, Type B2 | 100% discharged to outdoor space | | | | | | |
| Vertical laminar flow clean bench | 100% discharged to indoor space | Operation of non-hazardous bacteria, PIVAS | Positive pressure (Supply air to space outside of cabinet) | High efficiency | x | ✓ | x |
| Horizontal laminar flow clean bench | 100% discharged to indoor space | | | | | | |


Model Selection Guide for Biological Safety Cabinet

| | Applications | Class II, Type A2 | Class II, Type A2 + Discharge Ducting | Class II, Type B2 |
|---|---|-------------------|---------------------------------------|-------------------|
| Biotechnology | Sterilized culture medium preparation | ✓ | ✓ | ✓ |
| | Non-biohazard culture medium preparation | ✓ | ✓ | ✓ |
| | Culture | ✓ | ✓ | ✓ |
| | Non-biohazard tissue culture | ✓ | ✓ | ✓ |
| | Tissue culture | ✓ | ✓ | ✓ |
| | Plant tissue culture | ✓ | ✓ | ✓ |
| | Blood composition analysis | ✓ | ✓ | ✓ |
| | Human tissue research | ✓ | ✓ | ✓ |
| | PCR | ✓ | ✓ | ✓ |
| Microorganism | Sterilized culture medium preparation | ✓ | ✓ | ✓ |
| | Non-biohazard culture medium preparation | ✓ | ✓ | ✓ |
| | Culture | ✓ | ✓ | ✓ |
| | Odorless substance culture | | ✓ | ✓ |
| | Non-biohazard culture | ✓ | ✓ | ✓ |
| | Isolated clinical specimen | ✓ | ✓ | ✓ |
| | Blood analysis | ✓ | ✓ | ✓ |
| | QA/QC | ✓ | ✓ | ✓ |
| | Non-volatile toxic substance staining | ✓ | ✓ | ✓ |
| | Trace-volatile toxic substance staining | | ✓ | ✓ |
| | Non-volatile substance radioisotope labelling | ✓ | ✓ | ✓ |
| Trace-volatile substance radioisotope labelling | | ✓ | ✓ | |
| Medicine | Anticancer drug preparation | | ✓ | ✓ |
| | Trace-volatile substance preparation | | ✓ | ✓ |
| Routine Research | Cell/tissue fixation/staining | | ✓ | ✓ |
| | Toxic powder/suspended solids | ✓ | ✓ | ✓ |

NSF Series Biological Safety Cabinet

Scope of Application

Professional air purification equipment suitable for pharmaceuticals, medical and health industries, scientific research laboratories of universities and colleges as well as other related fields.



Product Advantages



Multiple Voltage Options, Suitable for Many Countries and Regions

Full voltage coverage (100-230V 50/60Hz), suitable for a wide range of countries and regions.



Microprocessor Control System

- Intuitive and informative interactive digital LCD display.
- The hot-ball anemometer monitors the downflow and inflow air speed of the safety cabinet in real time and compares it with the standard air speed. The rotating speed of the fan is adjusted through the microprocessor system to maintain the constant air speed of the safety cabinet.
- Real-time display of operational information and parameters including downflow air velocity, flow rate, temperature, humidity, positive pressure, negative pressure, fan cumulative running time and filter remaining service life.
- One button UV lamp timer function, allows users to set 0 to 24 hours of automatic on/off time.



Ultra-low Noise, Uniform Airflow

Designed with Dual DC centrifugal fans, combined with an innovative air distribution system, with lower noise, and more uniform air flow.



Superior Filter, Multiple Protection

- ULPA is made of moisture-proof and flame-retardant glass fiber filter paper which can intercept 99.9995% solid particles with a diameter of 0.12 μm to ensure high cleanliness of air supply flow and exhaust flow.
- Perfect air distribution design, no turbulence in the working area.
- Sound and light alarm function for abnormal parameters.



One-piece Welded Cabinet Structure, Leak Proof

Prevention of leakage performance of dangerous factors conforms to NSF specification.

NSF Series Biological Safety Cabinet

N Series Class II Type A2 Biological Safety Cabinet

Specifications

| Model | HR900-IIA2-N | HR1200-IIA2-N | HR1500-IIA2-N | HR1800-IIA2-N | | |
|-----------------------------------|--|--|---------------------------|---------------------------|--------------------|--------------------|
| Power Supply(V/Hz) | 100-230/50/60 | 100-230/50/60 | 100-230/50/60 | 100-230/50/60 | | |
| Power (VA) | 1300 | 1300 | 1400 | 1400 | | |
| Power of Blower (W) | DC 190W 112W | DC 120W 112W | DC 190W 112W | DC 120W 112W | | |
| Airflow Circulation | 70% Downflow, 30% Exhaust | 70% Downflow, 30% Exhaust | 70% Downflow, 30% Exhaust | 70% Downflow, 30% Exhaust | | |
| Main Filter Typical Efficiency | ULPA, U15,99.9995%@0.12um | ULPA, U15,99.9995%@0.12um | ULPA, U15,99.9995%@0.12um | ULPA, U15,99.9995%@0.12um | | |
| Exhaust Filter Typical Efficiency | ULPA, U15,99.9995%@0.12um | ULPA, U15,99.9995%@0.12um | ULPA, U15,99.9995%@0.12um | ULPA, U15,99.9995%@0.12um | | |
| Filter's Brand | AAF | AAF | AAF | AAF | | |
| Downflow Velocity (m/s) | 0.35 | 0.35 | 0.35 | 0.35 | | |
| Inflow Velocity (m/s) | 0.53 | 0.53 | 0.53 | 0.53 | | |
| Lighting Intensity (Lux) | 946 | 1323 | TBD | 1292 | | |
| Sound Level (dB(A)) | 67 | 64 | TBD | 67 | | |
| Dimensions | Net/Gross Weight | kg | 225/280 | 275/335 | 325/376 | 375/460 |
| | | lbs | 496/584 | 617/750 | 716.5/804.7 | 827/1014 |
| | Internal Dimensions (W*D*H) | mm | 930*600*650 | 1230*600*650 | 1530*600*650 | 1830*600*650 |
| | | in | 36.6*23.6*25.6 | 48.4*23.6*25.6 | 60.2*23.6*25.6 | 72.0*23.6*25.6 |
| | External Dimensions (W*D*H) (Width without Armrest) | mm | 1040*850(790)*2160 | 1340*850(790)*2160 | 1640*850(790)*2160 | 1940*850(790)*2160 |
| | | in | 40.9*33.5(31.1)*85 | 52.8*33.5(31.1)*85 | 64.6*33.5(31.1)*85 | 76.4*33.5(31.1)*85 |
| | Packing Dimensions (W*D*H) | mm | 1105*935*1720 | 1435*945*1700 | 1705*940*1720 | 2038*945*1700 |
| in | | 43.5*36.8*67.7 | 56.5*37.2*66.9 | 67.1*37.0*67.7 | 80.2*37.2*66.9 | |
| Support Stand (mm) | 680mm (standard), 680-900mm adjustable height (optional) | | | | | |
| Sash Opening (mm) | 200 (Max 460) | | | | | |
| Container Load (20'/40'/40'H) | 12/24/24 | | | | | |
| Alarms Functions | Alarm | Sound and Flash | | | | |
| | Closing Fan Alarm After Door Opening | Y | | | | |
| | Door Ajar | Y | | | | |
| | Abnormal Inflow | Y | | | | |
| | Abnormal Downflow | Y | | | | |
| | Door Open More than Limit | Y | | | | |
| | Filter Clogging | Y | | | | |
| | Filter Damage | Y | | | | |
| | Filter Lifetime Lack | Y | | | | |
| | UV Light Lifetime Lack | Y | | | | |
| Accessories | UV Lamp | Y | | | | |
| | Socket | 2 sockets (optional GFCI socket) | | | | |
| | Valve Port | Optional | | | | |
| | Side Wall Service Taps | Optional gas tap, vaccum tap, compressed air tap and water tap | | | | |
| | Antimicrobial Coating | Optional | | | | |
| | 316SS Work Surface | Optional | | | | |
| | Work Surface | One-piece | | | | |
| | Touch Screen | Optional | | | | |
| | Motorized Window | N | | | | |
| | Foot Switch | N | | | | |
| | Armrest | Y | | | | |
| | Exhaust Duct Connection Kits | Optional | | | | |
| | Others | Certification | UL, NSF | | | |

*Haier Biomedical reserves the right to change products and specifications without prior notice

Scope of Application

Designed to deliver airflow performance with exceptional precision for cell culture applications and other sensitive work requiring assured safety for product, personal and environment, especially for Bio-pharma and advanced research laboratories.

Dual DC fans - independent air supply and exhaust system to provide a uniform airflow and excellent energy efficiency, not affected by voltage fluctuations.

The ergonomic 10° inclination of the cabinet makes the operation more comfortable and the vision clearer

Equipped with an U15 ULPA: filter to provide 10x Filtration efficiency of HEPA filter, reaching the US standard Fed STD 209E class 1

Two Electrical Waterproof sockets with timing technology allows user to program timed on/off. Optional CFCI sockets

Intuitive Touch screen ① Digital display of running status, abnormal operation conditions alarm, filter and UV lamp service life warning function, UV lamp timing function and interlock technology

One-piece arm rest, comfortable to use, relieve fatigue.

Adjustable height supporter (680-900mm) (optional) ②

Universal casters

HR900-IIA2-N (NSF/ANSI 49 ONLY)
HR1200-IIA2-N
HR1500-IIA2-N
HR1800-IIA2-N

① Optional for NSF standard models

② Adjustable supporter is optional for NSF

N Series Class II Type A2 Biological Safety Cabinet

Product Advantages



Microprocessor Control System

The hot-ball anemometer monitors the downflow and inflow wind speed of the safety cabinet in real time and compares it with the standard wind speed. The rotating speed of the fan is adjusted through the microprocessor system to maintain the constant wind speed of the safety cabinet. Intuitive digital touch screen ①, real-time display of operational information and parameters including airflow velocity, temperature, humidity, positive pressure, negative pressure, fan cumulative running time and filter remaining service life. One button UV lamp timer function, allows users to set 0 to 24 hours of automatic on/off time.

① Optional for NSF standard models



Superior Filter, Multiple Protection

ULPA is made of moisture-proof and flame-retardant glass fiber filter paper which can intercept 99.9995% solid particles with a diameter of 0.12 μm to ensure high cleanliness of air supply flow and exhaust flow with 10x Filtration efficiency of HEPA filter.




Ultra-Low Noise, Uniform Airflow

Designed with Dual DC centrifugal fans, combined with an innovative air distribution system, with lower noise, and more uniform air flow.

Optional Components

Segmented Work Surfaces ①



Foot Switch ②



Exhaust Hood and Pipe



316 Stainless Steel Workbench and Liner




Valve Port



① CE12469 standard models only
② CE12469 standard models only


Product Features

Intuitive and User-friendly 10" Touchscreen Controller



- Cabinet parameters with understandable graphical illustration
- Real Time display of downflow and inflow velocity
- Scheduled UV appointment for most efficient decontamination
- Company common menu bar from Haier Biomedical
- Angled down for easy reach and viewing

Dual Wind Speed Sensors



Dual wind speed sensors-real-time monitoring the wind speed in the working area and exhaust area of the safety cabinet. Our microprocessor system will automatically adjust wind speed according to safety standard.

Wear-resistant and Easy to Clean



- The inner liner coating adopts stainless steel one-piece structure design, 12mm large arc transition, no cleaning dead angle.
- Operating table adopts high quality 304 stainless steel, without any fixed screws, will not form the accumulation of pollutants.
- The cabinet is made of cold rolled steel plate with corrosion resistant epoxy resin powder sprayed on the surface, which is resistant to acid and alkali corrosion.

Caster and Foot



4 universal casters +4 feet design, easy to move and lock leveling.

Platform Arm Rest



Platform arm rest frame, comfortable to use, relieve fatigue.

N Series Class II Type A2 Biological Safety Cabinet

Specifications

| | | NSF/ANSI 49 | | | | EN12469 | | |
|---|--|--|--|--|--|--|--|--|
| General Specifications | Working Voltage&Frequency (V/Hz) | 100-230/50/60 | 100-230/50/60 | 100-230/50/60 | 100-230/50/60 | 220-240/50/60 | 220-240/50/60 | 220-240/50/60 |
| | Nominal Power Consumption (W) | 200 | 200 | 240 | 440 | 200 | 240 | 440 |
| | Heat Load (BTU/Hr) | 682 | 682 | 818 | 1500 | 682 | 818 | 1500 |
| | Power of Blower (W) | DC 120, DC 112 | DC 120, DC 112 | DC 190, DC 112 | DC 120, DC 112 | DC 120, DC 112 | DC 190, DC 112 | DC 120, DC 112 |
| | Airflow Circulation | 70% Downflow,30% Exhaust | 70% Downflow,30% Exhaust | 70% Downflow,30% Exhaust | 70% Downflow,30% Exhaust | 70% Downflow,30% Exhaust | 70% Downflow,30% Exhaust | 70% Downflow,30% Exhaust |
| | Main Filter Typical Efficiency | ULPA,U15.99.9995%@0.12um | ULPA,U15.99.9995%@0.12um | ULPA,U15.99.9995%@0.12um | ULPA,U15.99.9995%@0.12um | ULPA,U15.99.9995%@0.12um | ULPA,U15.99.9995%@0.12um | ULPA,U15.99.9995%@0.12um |
| | Exhaust Filter Typical Efficiency | ULPA,U15.99.9995%@0.12um | ULPA,U15.99.9995%@0.12um | ULPA,U15.99.9995%@0.12um | ULPA,U15.99.9995%@0.12um | ULPA,U15.99.9995%@0.12um | ULPA,U15.99.9995%@0.12um | ULPA,U15.99.9995%@0.12um |
| | Filter's Brand | AAF | AAF | AAF | AAF | AAF | AAF | AAF |
| | Downflow Velocity (m/s) | 0.35 | 0.35 | 0.35 | 0.35 | 0.30 | 0.30 | 0.30 |
| | Inflow Velocity (m/s) | 0.53 | 0.53 | 0.53 | 0.53 | 0.45 | 0.45 | 0.45 |
| | Inflow Volume (m³/h) | 390 | 516 | 642 | 768 | 438 | 545 | 652 |
| | Lighting Intensity (Lux) | 946 | 1323 | 1248 | 1292 | 1233 | 1313 | 1212 |
| | Fan | DC+DC | DC+DC | DC+DC | DC+DC | DC+DC | DC+DC | DC+DC |
| | Fan's Brand | EBM | EBM | EBM | EBM | EBM | EBM | EBM |
| Solid or Glass Sides | Solid | Solid | Solid | Solid | Solid | Solid | Solid | |
| Sash Opening (mm) | 200(Max 460) | 200(Max 460) | 200(Max 460) | 200(Max 460) | 200(Max 460) | 200(Max 460) | 200(Max 460) | |
| Sound Level(dB(A)) EU units measured according to EN 12469, NSF units measured according to NSF/ANSI 49 | 67 | 64 | 63 | 67 | 58.9 | 60.5 | 62.8 | |
| Dimensions | Net/Gross Weight | 225/280 496/617 | 275/335 617/750 | 325/376 716.5/804.7 | 375/460 827/1014 | 275/335 617/750 | 325/376 716.5/804.7 | 375/460 827/1014 |
| | Internal Dimensions (W*D*H) | 930*600*650 36.6*23.6*25.6 | 1230*600*650 48.4*23.6*25.6 | 1530*600*650 60.2*23.6*25.6 | 1830*600*650 72.0*23.6*25.6 | 1230*600*650 48.4*23.6*25.6 | 1530*600*650 60.2*23.6*25.6 | 1830*600*650 72.0*23.6*25.6 |
| | External Dimensions (W*D*H) (Width without Arm Rest) | 1040*850(790)*2160 40.9*33.5(31.1)*85 | 1340*850(790)*2160 52.8*33.5(31.1)*85 | 1640*850(790)*2160 64.6*33.5(31.1)*85 | 1940*850(790)*2160 76.4*33.5(31.1)*85 | 1340*850(790)*2160 52.8*33.5(31.1)*85 | 1640*850(790)*2160 64.6*33.5(31.1)*85 | 1940*850(790)*2160 76.4*33.5(31.1)*85 |
| | Packing Dimensions (W*D*H) | 1105*935*1720 43.5*36.8*67.7 | 1435*945*1700 56.5*37.2*66.9 | 1705*940*1720 67.1*37.0*67.7 | 2038*945*1700 80.2*37.2*66.9 | 1435*945*1700 56.5*37.2*66.9 | 1705*940*1720 67.1*37.0*67.7 | 2038*945*1700 80.2*37.2*66.9 |
| | Support Stand (mm) | 680mm(standard), 680-900mm adjustable height(optional) | 680mm(standard), 680-900mm adjustable height(optional) | 680mm(standard), 680-900mm adjustable height(optional) | 680mm(standard), 680-900mm adjustable height(optional) | 680-900mm adjustable height (standard) | 680-900mm adjustable height (standard) | 680-900mm adjustable height (standard) |
| | Container Load (20'/40'/40'H) | 8/16/16 | 8/16/16 | 6/12/12 | 6/12/12 | 8/16/16 | 6/12/12 | 6/12/12 |
| Alarm Functions | Alarm | Sound and Flash | Sound and Flash | Sound and Flash | Sound and Flash | Sound and Flash | Sound and Flash | Sound and Flash |
| | Closing Fan Alarm After Door Opening | Y | Y | Y | Y | Y | Y | Y |
| | Door Ajar | Y | Y | Y | Y | Y | Y | Y |
| | Abnormal Inflow | Y | Y | Y | Y | Y | Y | Y |
| | Abnormal Downflow | Y | Y | Y | Y | Y | Y | Y |
| | Door Open More than Limit | Y | Y | Y | Y | Y | Y | Y |
| | Filter Clogging | Y | Y | Y | Y | Y | Y | Y |
| | Filter Damage | Y | Y | Y | Y | Y | Y | Y |
| | Filter Lifetime Lack | Y | Y | Y | Y | Y | Y | Y |
| | UV Light Lifetime Lack | Y | Y | Y | Y | Y | Y | Y |
| Accessories | Service Ports (Port Holes) | Y | Y | Y | Y | Y | Y | Y |
| | UV Lamp | Y | Y | Y | Y | Y | Y | Y |
| | Socket | 2 Sockets (optional GFCI socket) | 2 Sockets (optional GFCI socket) | 2 Sockets (optional GFCI socket) | 2 Sockets (optional GFCI socket) | 2 Sockets | 2 Sockets | 2 Sockets |
| | Side Wall Service Taps | Optional gas tap, vaccum tap, compressed air tap and water tap | Optional gas tap, vaccum tap, compressed air tap and water tap | Optional gas tap, vaccum tap, compressed air tap and water tap | Optional gas tap, vaccum tap, compressed air tap and water tap | Optional gas tap, vaccum tap, compressed air tap and water tap | Optional gas tap, vaccum tap, compressed air tap and water tap | Optional gas tap, vaccum tap, compressed air tap and water tap |
| | Antimicrobial Powder Coating | Optional | Optional | Optional | Optional | Optional | Optional | Optional |
| | 316SS Work Surface | Optional | Optional | Optional | Optional | Optional | Optional | Optional |
| | Work Surface | One-piece | One-piece | One-piece | One-piece | One-piece | One-piece | One-piece |
| | Touch Screen | Optional | Optional | Optional | Optional | Y | Y | Y |
| | Motorized Window | N | N | N | N | Optional | Optional | Optional |
| | Foot Switch | N | N | N | N | Optional | Optional | Optional |
| Others | Armrest | Y | Y | Y | Y | Y | Y | Y |
| | Exhaust Duct Connection Kits | Optional | Optional | Optional | Optional | Optional | Optional | Optional |
| Certification | UL NSF | UL NSF | UL NSF | UL NSF | CE/EN12469/TUVMark | CE/EN12469/TUVMark | CE/EN12469/TUVMark | |

Biological Safety Cabinet (Single/Double Exhaust Filter)

Scope of Application

Suitable for microbiology, biomedicine, biosafety laboratories and other laboratories. It offers three levels of protection-operator, product and environment



Product Advantages



Ultra-efficient Filtration

High efficiency moisture proof and flame-retarded filter achieving a retention efficiency of 99.9995% for 0.12 micron solid particle system



Intelligent Constant Air Velocity System

Designed with intelligent constant air speed technology, which can ensure that the working area could have its downflow and inflow velocities in line with the standard requirements



Side Glass Windows

Tempered glass side walls increase visibility and prevents the operator from experiencing a "boxed-in" feeling



Multiple Alarm Functions

Audible and visual alarm functions: filter and UV end-of-life alerts, fan turned-off after door opening alert and door open alarm



Ergonomic Design

- A comfortable platform-type armrest can reduce hand and arm fatigue
- 10-degree angle inclination design, in accordance with ergonomic principles for more comfortable operation
- With height adjustable stand, users can adjust the height of the work bench to requirements



Easy to Clean

The armrest can be removed and the unique drop-down front sash can be removed in seconds to enable easy cleaning and full access to the cabinet interior work space for thorough cleaning/disinfection



Lower Noise and Airflow Uniformity

DC&EC fans operate with lower noise and deliver excellent air flow uniformity



Patented IP-44 Rated Socket with Timer

The sockets can be programmed to supply power at a specified time to meet the users' experiment schedules



UV Lamp

Patented one-button UV lamp technology allows users to automatically activate or deactivate based on usage habits or at specific sterilization intervals to minimise waiting time



PIR Detection

In intelligent mode, the PIR sensor will detect when the operating area has been clear of personnel for more than 15 minutes and automatically switch to LNS green-saving mode to reduce the noise level, save energy and prolong the service life of filter

Biological Safety Cabinet (Single/Double Exhaust Filter)

Specifications

| Model | | HR900-IIA2-S | HR1200-IIA2-S | HR900-IIA2-D | HR1200-IIA2-D | |
|-------------------------------|---|--|--|--|--|----------------------|
| General Specification | Power Supply(V/Hz) | 220-240/50/60 | 220-240/50/60 | 220-240/50/60 | 220-240/50/60 | |
| | Power (VA) | 1500 | 1600 | 1500 | 1600 | |
| | Power of Blower (W) | DC 120W EC170W | DC 120, DC 112 | DC 120W EC170W | DC 190, DC 170 | |
| | Airflow Circulation | 70% Downflow, 30% Exhaust | 70% Downflow, 30% Exhaust | 70% Downflow, 30% Exhaust | 70% Downflow, 30% Exhaust | |
| | Main Filter Typical Efficiency | ULPA, U15, 99.9995% @ 0.12um | HEPA, H14, 99.995% @ 0.3um | ULPA, U15, 99.9995% @ 0.12um | ULPA, U15, 99.9995% @ 0.12um | |
| | Exhaust Filter Typical Efficiency | ULPA, U15, 99.9995% @ 0.12um | HEPA, H14, 99.995% @ 0.3um | ULPA x2, U15, 99.9995% @ 0.12um | HEPA x2, H14, 99.995% @ 0.3um | |
| | Exhaust Filter | Single | Single | Dual | Dual | |
| | Filter's Brand | AAF | AAF | AAF | AAF | |
| | Downflow Velocity (m/s) | 0.30 | 0.30 | 0.30 | 0.30 | |
| | Inflow Velocity (m/s) | 0.45 | 0.45 | 0.45 | 0.45 | |
| Lighting Intensity (Lux) | ≥1300 | ≥1000 | ≥1300 | ≥1000 | | |
| Sound Level (dB(A)) | 57.1 | 60 | 59.7 | 60 | | |
| Dimensions | Net/Gross Weight | kg | 290/320 | 320/339 | 290/320 | 320/339 |
| | | lbs | 639/705 | 705.5/747.4 | 639/705 | 705.5/747.4 |
| | Internal Dimensions (W*D*H) | mm | 936*620*635 | 1310*620*635 | 936*620*635 | 1310*620*635 |
| | | in | 36.9*24.4*25 | 51.6*24.4*25 | 36.9*24.4*25 | 51.6*24.4*25 |
| | External Dimensions (W*D*H) (Width without Armrest) | mm | 1002*856(796)*1485 | 1380*856(796)*2164 | 1002*856(796)*1485 | 1380*856(796)*1485 |
| | | in | 39.5*33.7(31.3)*58.5 | 54.3*33.7(31.3)*85.2 | 39.5*33.7(31.3)*58.5 | 54.3*33.7(31.3)*58.5 |
| | Packing Dimensions (W*D*H) | mm | 1185*925*1720 | 1475*935*1720 | 1105*932*1730 | 1475*935*1720 |
| in | | 46.7*36.4*67.7 | 58.1*36.8*67.7 | 43.5*36.7*68.1 | 58.1*36.8*67.7 | |
| Support Stand (mm) | No, 680-900mm adjustable height (optional) | | 680-900mm adjustable height (standard) | No, 680-900mm adjustable height (optional) | No, 680-900mm adjustable height (optional) | |
| Sash Opening (mm) | 200 (Max 490) | | 200 (Max 490) | 200 (Max 490) | 200 (Max 490) | |
| Container Load (20"/40"/40"H) | 11/21/21 | | 8/16/16 | 11/21/21 | 8/16/16 | |
| Alarms Functions | Alarm | Sound and Flash | Sound and Flash | Sound and Flash | Sound and Flash | |
| | Closing Fan Alarm After Door Opening | Y | Y | Y | Y | |
| | Door Ajar Alarm | Y | Y | Y | Y | |
| | Abnormal Inflow Alarm | Y | Y | Y | Y | |
| | Abnormal Downflow Alarm | Y | Y | Y | Y | |
| | Door Open More than Limit Alarm | Y | Y | Y | Y | |
| | Filter Clogging Alarm | Y | Y | Y | Y | |
| | Filter Damage Alarm | Y | Y | Y | Y | |
| | Filter Lifetime Lack Alarm | Y | Y | Y | Y | |
| | UV Light Lifetime Lack Alarm | Y | Y | Y | Y | |
| Accessories | UV Lamp | Y | Y | Y | Y | |
| | Socket | 2 sockets | 2 sockets | 2 sockets | 2 sockets | |
| | Valve Port | Optional | Optional | Optional | Optional | |
| | Side Wall Service Taps | Optional gas tap, vacuum tap, compressed air tap and water tap | Optional gas tap, vacuum tap, compressed air tap and water tap | Optional gas tap, vacuum tap, compressed air tap and water tap | Optional gas tap, vacuum tap, compressed air tap and water tap | |
| | Antimicrobial Coating | Optional | Optional | Optional | Optional | |
| | 316SS Work Surface | Optional | Optional | Optional | Optional | |
| | Work Surface | Three-piece | Four-piece | Three-piece | Three-piece | |
| | Touch Screen | N | N | N | N | |
| | Motorized Window | Optional | Optional | Optional | Optional | |
| | Foot Switch | Optional | Optional | Optional | Optional | |
| Others | Armrest | Y | Y | Y | Y | |
| | Exhaust Duct Connection Kits | Optional | Optional | Optional | Optional | |
| Certification | CE, UKCA | CE, UKCA, TUV | CE, UKCA | CE, UKCA, TUV | | |

*Haier Biomedical reserves the right to change products and specifications without prior notice

Biological Safety Cabinet

Scope of Application

This Class II microbiological safety cabinet is designed to protect the operator, laboratory environment and samples from being exposed to the infective aerosol produced from samples with bacteria strains, diagnostic materials, and other infective substances. It provides the operator with comfortable and safer working conditions. It is widely used in medical health, disease prevention, food safety, biological pharmacy and environment monitoring.

HR1200-IIA2

- ULPA filter (U15)
- Forward centrifugal fan
- Intelligent constant airflow velocity
- Patented airflow technology Interlocked internal & external fans Airflow valve interlocking control
- Stainless steel armrest
- Universal caster with built-in threaded supporting leg
- Digital display of operation parameters Filter end-of-life (EOL) reminder UV lamp end-of-life (EOL) reminder UV lamp appointment, patented technology UV lamp interlocking, patented technology
- Abnormal operation condition alarm Key component failure alarm
- Cabinet body with 10° inclination
- PIR detection
- Drop-down front glass window
- Patented IP-44 rated socket with timer Power supply plug with leakage protection
- V-shaped air inlet

Biological Safety Cabinet

Product Advantages



V-shaped Air Inlet

The V-shaped air inlet can prevent the samples or arms of operator from blocking the air flow.

The work surface can be easily lifted using the handles for cleaning purposes



One-touch UV Lamp Operation

UV lamp records and remembers users' setting and habits and can be preset with a startup delay with one-key operation for ease and convenience



Energy Conservation

- PIR detection

In intelligent mode, the PIR sensor will detect when the operating area has been clear of personnel for more than 15 minutes and automatically switch to LNS green-saving mode to reduce the noise level, save energy and prolong the service life of filter



Ergonomic

- 10° inclination design of cabinet body

The front operation interface has an ergonomic design of 10 ° inclination for ensuring more comfortable operation

- Stainless steel armrest

A comfortable platform-type armrest can reduce hand and arm fatigue

Easy to Clean

The armrest can be removed and the unique drop-down front sash can be removed in seconds to enable easy cleaning and full access to the cabinet interior work space for thorough cleaning/disinfection



Intelligent

Constant airflow velocity

The hot-bulb airflow velocity transducer performs real-time monitoring of the air velocity of the working area, compares it with the standard air velocity, and keeps a constant air velocity in the cabinet by regulating the fan speed with a microcomputer system



IP 44 Rated Power Sockets with Timer

The sockets can be programmed to supply power at a specified time to meet the users' demands for timing of experiments



Safe

- Abnormal operation condition alarm

Audible and visual alarms in form of voice or text will be present when air turbulence level exceeds 20% and door height (high or low) or work area temperature exceeds limits

- Patented technology: filter end-of-life reminder

Pressure transducer monitors the resistance variation of filter to determine the remaining life of filter and will remind the user by warning when the remaining life is below 10%

- Patented technology: UV lamp end-of-life reminder

The microcomputer will add up the service time of UV lamp, and will remind by warning the user to replace the UV lamp when its remaining life is less than 10%




Offers a Broad Range of Sizes to Fit any Available Laboratory Space

Mini Biological Safety Cabinet

Scope of Application

Small size microbiological safety cabinet for single person operations and small working spaces, suitable for testing stations within warehouse, airports and mobile laboratories.



HR700-IIA2

Mini Biological Safety Cabinet

Product Advantages

Safe and Reliable
Adopts AAF filters, class 1 cleanliness to provide maximum protection for people, environment and samples

Interlocking Function
The product features an interlocking function between UV disinfection, fluorescent lamp, front sash and fan

LCD Display
LCD screen displays various parameters and service life of accessories in real time, and the operational condition of the equipment is clear at a glance

Alarms
Audible and visual alarm functions: hardware failure alarm, operation parameter overrun alarm, filter and UV lamp life warning

One-click Operation
UV lamp one-click reservation allows users to set 0 min to 24 hours of automatic on/off time and sterilization interval, reducing the waiting time

Uniform Airflow
The EBM fan, combined with Haier Biomedical's professional air distribution design, provides lower noise operation and a uniform air flow

Two Waterproof Sockets
Sockets include timing technology to allow users to program socket on/off times

Safety Certified
Complies with EN12469



Specifications


| Model | HR700-IIA2 | HR900-IIA2 | HR1200-IIA2 | HR1500-IIA2 |
|--------------------------------------|--|--|--|--|
| Power Supply(V/Hz) | 220-240/50 | 220-240/50 | 220-240/50 | 220-240/50 |
| Power (VA) | 1200 | 1500 | 1500 | 1900 |
| Power of Blower (W) | AC 210 | AC L=330,M=465,H=735 | AC L=330,M=465,H=735 | AC 650 |
| Airflow Circulation | 70% Downflow,30% Exhaust | 70% Downflow,30% Exhaust | 70% Downflow,30% Exhaust | 70% Downflow,30% Exhaust |
| Main Filter Typical Efficiency | HEPA,H14,99.995%@0.3um | ULPA,U15,99.9995%@0.12um | ULPA,U15,99.9995%@0.12um | ULPA,U15,99.9995%@0.12um |
| Exhaust Filter Typical Efficiency | ULPA,U15,99.9995%@0.12um | ULPA,U15,99.9995%@0.12um | HEPA,H14,99.995%@0.3um | HEPA,H14,99.995%@0.3um |
| Exhaust Filter | Single | Single | Single | Single |
| Filter's Brand | AAF | AAF | AAF | AAF |
| Downflow Velocity (m/s) | 0.30 | 0.33 | 0.34 | 0.31 |
| Inflow Velocity (m/s) | 0.53 | 0.55 | 0.55 | 0.55 |
| Lighting Intensity (Lux) | ≥780 | ≥900 | ≥900 | ≥900 |
| Sound Level (dB(A)) | <65 | 60 | 60 | 62 |
| Net/Gross Weight | kg | 100/130 | 290/310 | 320/339 |
| | lbs | 220.4/286.52 | 639.3/683.4 | 705.5/747.4 |
| Internal Dimensions (W*D*H) | mm | 600*550*540 | 920*620*650 | 1220*620*650 |
| | in | 23.62*21.65*21.26 | 36.2*24.4*25.6 | 48.0*24.4*25.6 |
| External Dimensions (W*D*H) | mm | 700*720*1200 | 1080*845(790)*2160 | 1380*845(790)*2160 |
| | in | 27.56*28.35*47.24 | 42.5*33.3(31.1)*85.0 | 54.3*33.3(31.1)*85.0 |
| Packing Dimensions (W*D*H) | mm | 800*810*1385 (without stand) 825*820*1495 (with stand) | 1185*925*1720 | 1475*935*1720 |
| | in | 31.5*31.9*54.5 (without stand) 32.5*32.3*58.9 (with stand) | 46.7*36.4*67.7 | 58.1*36.8*67.7 |
| Support Stand (mm) | No,700mm (optional) | 680-900mm adjustable height (standard) | 680-900mm adjustable height (standard) | 680-900mm adjustable height (standard) |
| Sash Opening (mm) | 200 (Max) | 200 (Max 490) | 200 (Max 490) | 200 (Max 490) |
| Container Load (20"/40"/40"H) | 14/28/28 | 11/21/21 | 8/16/16 | 6/12/12 |
| Alarm | Sound and Flash | Sound and Flash | Sound and Flash | Sound and Flash |
| Closing Fan Alarm After Door Opening | Y | Y | Y | Y |
| Door Ajar Alarm | Y | Y | Y | Y |
| Abnormal Inflow Alarm | Y | Y | Y | Y |
| Abnormal Downflow Alarm | Y | Y | Y | Y |
| Door Open More than Limit Alarm | Y | Y | Y | Y |
| Filter Clogging Alarm | Y | Y | Y | Y |
| Filter Damage Alarm | Y | Y | Y | Y |
| Filter Lifetime Lack Alarm | Y | Y | Y | Y |
| UV Light Lifetime Lack Alarm | Y | Y | Y | Y |
| UV Lamp | Y | Y | Y | Y |
| Socket | 2 sockets | 2 sockets | 2 sockets | 2 sockets |
| Valve Port | Optional | Standard with 2 valve port | Standard with 2 valve port | Standard with 2 valve port |
| Side Wall Service Taps | Optional gas tap, vaccum tap, compressed air tap and water tap | Standard gas tap and vaccum tap; optional compressed air tap and water tap | Standard gas tap and vaccum tap; optional compressed air tap and water tap | Standard gas tap and vaccum tap; optional compressed air tap and water tap |
| Antimicrobial Coating | Optional | Optional | Optional | Optional |
| 316SS Work Surface | Optional | Optional | Optional | Optional |
| Work Surface | One-piece | One-piece | One-piece | One-piece |
| Touch Screen | N | N | N | N |
| Electric Front Sash | Y | Optional | Optional | Optional |
| Foot Switch | Optional | Optional | Optional | Optional |
| Armrest | N | Y | Y | Y |
| Exhaust Duct Connection Kits | N | Optional | Optional | Optional |
| Others | Certification | CE, TUV | CE, UKCA, TUV | CE, UKCA, TUV |

*Haier Biomedical reserves the right to change products and specifications without prior notice

Classic Series Biological Safety Cabinet (Type A2)

Scope of Application

Professional partial air purification equipment, it is suitable for cell biology, microbiology, biomedicine, biosafety and other related laboratories.



Product Advantages

- Patented air flow blocking technology is designed at the upper edge and both edges of front window to eliminate the exposure of microorganisms
- The internal walls on three sides of operation area is constructed by a single plate, and the 12mm arc angle corner for optimal cleaning
- Quality 304 stainless steel work surface without screws, no accumulation of contaminant
- Digital LCD screen
- Dismountable air in-flow plate, easy to clean and sterilize
- UV lamp one-click reservation allows users to set 0 min to 24 hours of automatic on/off time and sterilization interval, reducing the waiting time
- Professional**
Digital display of operating parameters
Real-time digital display of down flow, inflow, exhaust volume, filter remaining life, UV lamp remaining life, negative pressure and positive pressure
- Real-time display of key parameters: down-flow velocity, inflow velocity, airflow volume, static pressure, negative pressure, accumulative running time off an accumulative running time of UV lamp, and remaining life of filter
- The volume of liquid tank is over 4L, designed with outlet valve for convenient cleaning and maintenance
- Audible and visual alarms for abnormal parameters
- Safe**
 - Abnormal operation condition alarm**
Audible and visual alarms in form of voice or text will be present when air turbulence level exceeds 20% and door height (high or low)
 - Filter end-of-life reminder**
Pressure transducer monitors the resistance variation of filter to determine the remaining life of filter and reminds the user by warning when the remaining life is below 10%
 - UV lamp end-of-life reminder**
The microcomputer records service times of the UV lamp and will alert the user to replace the UV lamp when its remaining life is less than 10%
 - Interlocking feature to ensure high safety and reliability**
Patented technology: UV lamp interlocking control UV lamp, sash, fan and interior light operation are interlocked together to protect against harmful UV rays and to prevent leakage of microorganisms

Classic Series Biological Safety Cabinet (Type A2)

Specifications

| Model | HR30-IIA2 | HR40-IIA2 | HR40-IIA2 | | |
|-------------------------------|--|--|--|--|----------------|
| General Specification | Power Supply(W/Hz) | 220-240/50/60 | 115/60 | 220-240/50/60 | |
| | Power (VA) | 1300 | 1300 | 1300 | |
| | Power of Blower (W) | AC 540/625 | AC 540/625 | AC 540/625 | |
| | Airflow Circulation | 70% Downflow,30% Exhaust | 70% Downflow,30% Exhaust | 70% Downflow,30% Exhaust | |
| | Main Filter Typical Efficiency | ULPA ,U15.99.9995%@0.12um | ULPA ,U15.99.9995%@0.12um | ULPA ,U15.99.9995%@0.12um | |
| | Exhaust Filter Typical Efficiency | HEPA ,H14.99.9995%@0.3um | HEPA ,H14.99.9995%@0.3um | HEPA ,H14.99.9995%@0.3um | |
| | Exhaust Filter | Single | Single | Single | |
| | Filter's Brand | AirePlus | AAF | AAF | |
| | Downflow Velocity (m/s) | 0.3 | 0.28 | 0.3 | |
| | Inflow Velocity (m/s) | 0.53 | 0.55 | 0.53 | |
| | Lighting Intensity (Lux) | ≥1100 | ≥1200 | ≥1200 | |
| | Sound Level (dB(A)) | 61 | 64 | 64 | |
| | Dimensions | Bench or Stand Only | Stand only | Stand only | Stand only |
| Net/Gross Weight | | kg | 220/248 | 258/305 | 258/305 |
| | | lbs | 485/546.7 | 568.8/672.4 | 568.8/672.4 |
| | | mm | | | |
| Internal Dimensions (W*D*H) | | in | 900*610*680 | 1167*610*680 | 1167*610*680 |
| | | mm | 35.4*24.0*26.8 | 45.9*24.0*26.8 | 45.9*24.0*26.8 |
| External Dimensions (W*D*H) | | in | 1100*820*2200 | 1360*820*2200 | 1360*820*2200 |
| | | mm | 43.3*32.3*86.6 | 53.5*32.3*86.6 | 53.5*32.3*86.6 |
| Packing Dimensions (W*D*H) | | in | 1185*945*1750 | 1445*945*1720 | 1445*945*1720 |
| | | mm | 46.7*37.2*68.9 | 56.9*37.2*67.7 | 56.9*37.2*67.7 |
| Support Stand (mm) | 680mm (standard), 680-900mm adjustable height (optional) | 680mm (standard), 680-900mm adjustable height (optional) | 680mm (standard), 680-900mm adjustable height (optional) | | |
| Sash Opening (mm) | 190 (Max 490) | 190 (Max 490) | 190 (Max 490) | | |
| Container Load (20"/40"/40"H) | 10/20/20 | 8/16/16 | 8/16/16 | | |
| Alarms Functions | Alarm | Sound and Flash | Sound and Flash | Sound and Flash | |
| | Closing Fan Alarm After Door Opening | Y | Y | Y | |
| | Door Ajar | Y | Y | Y | |
| | Abnormal Inflow | Y | Y | Y | |
| | Abnormal Downflow | Y | Y | Y | |
| | Door Open More than Limit | Y | Y | Y | |
| | Filter Clogging | Y | Y | Y | |
| | Filter Damage | Y | Y | Y | |
| | Filter Lifetime Lack | Y | Y | Y | |
| | UV Light Lifetime Lack | Y | Y | Y | |
| Accessories | UV Lamp | Y | Y | Y | |
| | Socket | 1 Socket | 1 Socket | 1 Socket | |
| | Valve Port | Standard with 2 Valve Port | Standard with 2 Valve Port | Standard with 2 Valve Port | |
| | Side Wall Service Taps | Standard gas tap and vaccum Tap; optional compressed air tap and water tap | Standard gas tap and vaccum Tap; optional compressed air tap and water tap | Standard gas tap and vaccum Tap; optional compressed air tap and water tap | |
| | Antimicrobial Coating | Optional | Optional | Optional | |
| | 316SS Work Surface | Optional | Optional | Optional | |
| | Work Surface | One-piece | One-piece | One-piece | |
| | Touch Screen | N | N | N | |
| | Motorized Window | N | N | N | |
| | Foot Switch | N | N | N | |
| Others | Armrest | Y | Y | Y | |
| | Exhaust Duct Connection Kits | Optional | Optional | Optional | |
| Certification | CE, TUV | / | CE, TUV | | |

*Haier Biomedical reserves the right to change products and specifications without prior notice

Classic Series Biological Safety Cabinet (Type B2)

Scope of Application

Professional partial air purification equipment, it is suitable for cell biology, microbiology, biomedicine, biosafety and other related laboratories.

Product Advantages



100% exterior exhaust



4m corrosion-resistant corrugated hose (standard)



Haier Biomedical exterior exhaust fan (optional)



Interlocking between safety cabinet and exterior exhauster, enables remote control of exterior exhauster parameters with safety cabinet



HR40-IIB2

Classic Series Biological Safety Cabinet (Type B2)

Animal Containment Workstation

Specifications

| Model | | HR40-IIB2 | |
|-------------------------------|--|--|----------------|
| General Specification | Power Supply(V/Hz) | 220-240/50/60 | |
| | Power (VA) | 1700 | |
| | Power of Blower (W) | AC 115 | |
| | Airflow Circulation | 100% Exhaust | |
| | Main Filter Typical Efficiency | ULPA ,U15.99.9995%@0.12um | |
| | Exhaust Filter Typical Efficiency | HEPA ,H14.99.995%@0.3um | |
| | Exhaust Filter | Single | |
| | Filter's Brand | AAF | |
| | Downflow Velocity (m/s) | 0.28 | |
| | Inflow Velocity (m/s) | 0.55 | |
| | Lighting Intensity (Lux) | ≥1200 | |
| | Sound Level (dB(A)) | 64 | |
| | Bench or Stand Only | Stand only | |
| Dimensions | Net/Gross Weight | kg | 252/308 |
| | | lbs | 555.6/679.0 |
| | Internal Dimensions (W*D*H) | mm | 1167*610*680 |
| | | in | 45.9*24.0*26.8 |
| | External Dimensions (W*D*H) | mm | 1360*820*2330 |
| | | in | 53.5*32.3*91.7 |
| | Packing Dimensions (W*D*H) | mm | 1445*945*1930 |
| | | in | 56.9*37.2*76.0 |
| Support Stand (mm) | 680mm (standard), 680-900mm adjustable height (optional) | | |
| Sash Opening (mm) | 190 (Max 490) | | |
| Container Load (20'/40'/40'H) | 8/16/16 | | |
| Alarms Functions | Alarm | Sound and Flash | |
| | Closing Fan Alarm After Door Opening | Y | |
| | Door Ajar | Y | |
| | Abnormal Inflow | Y | |
| | Abnormal Downflow | Y | |
| | Door Open More than Limit | Y | |
| | Filter Clogging | Y | |
| | Filter Damage | Y | |
| | Filter Lifetime Lack | Y | |
| | UV Light Lifetime Lack | Y | |
| Accessories | UV Lamp | Y | |
| | Socket | 1 Socket | |
| | Valve Port | Standard with 2 Valve Port | |
| | Side Wall Service Taps | Standard gas tap and vaccum tap; optional compressed air tap and water tap | |
| | Antimicrobial Coating | Optional | |
| | 316SS Work Surface | Optional | |
| | Work Surface | One-piece | |
| | Touch Screen | N | |
| | Motorized Window | N | |
| | Foot Switch | N | |
| Others | Armrest | Y | |
| | Exhaust Duct Connection Kits | N | |
| Certification | CE, TUV | | |

*Haier Biomedical reserves the right to change products and specifications without prior notice

Scope of Application

Applicable for laboratory animal operations (surgery, weighing, etc.) in experimental animal centers of scientific research institutions, hospitals, colleges and universities, or experimental animal research and development projects and breeding institutions.



HR1500-IIA2-DW

Product Advantages



Low Maintenance Costs

Experimental animal hair and debris are captured to eliminate the impact of hair debris and other debris on the service life of the fans and filters which ensures the cabinet is more durable



Sectioned Pre-filter

The interior of the sectioned pre-filter is easy to clean, and the pre-filter with a rear-mounted work surface is designed in sections, meaning the filter is easy to disassembled and sterilized as needed



The Dual System Ensures Maximum Safety

Each independent-controlled down flow and inlet airflow system is designed with dual fans, which ensures the stable airflow in the cabinet even if the pre-filter is 50% blocked. In the case of either damage to either fan, the air supply and exhaust system ensure safe air and protection



Stable Laminar Flow Offers Safety

The intelligent constant air speed design means that if abnormal air velocity is detected the fan adjust automatically to ensure the air speed stay within normal range

Product Features



Safety Interlock

Ultraviolet lamp, sash window, fans, and lamps are interlocked to eliminate hidden dangers of ultraviolet and microorganism leakage



Traceable

The work area in the cabinet is monitored by camera and the images are stored for future reference (optional)



Real-time Status Display

The digital screen displays in real time the flow rate of downflow and inflow airflow, the volume of discharged air, the temperature and humidity of the working area, the remaining service life of the filter, the UV lamp, the pressure in the negative pressure area, and positive pressure area. Users are alerted if abnormal operation is detected



Electric Foot Switch

When it is inconvenient to open the safety cabinet door with both hands, the electric front sash can be opened by foot switch (optional)

Laminar Flow Cabinet

Specifications

| Model | HR1500-IIA2-DW |
|---------------------------------|----------------|
| Power (W) | 1400 |
| Downflow Velocity (m/s) | 0.3 |
| Inflow Velocity (m/s) | 0.53 |
| Lighting Intensity (Lux) | ≥1100 |
| Net/Gross Weight (kg) | 350/400 |
| Interior Dimensions (W*D*H)(mm) | 1530*600*653 |
| Exterior Dimensions (W*D*H)(mm) | 1636*790*2170 |
| Packing Dimensions (W*D*H)(mm) | 1700*925*1715 |
| Container Load (20'/40'/40'H) | 6/12/12 |

*Haier Biomedical reserves the right to change products and specifications without prior notice

Scope of Application






A professional localized air purification equipment suitable for pharmaceuticals, medical and health, scientific research laboratories of universities and colleges, photoelectric/microelectronics manufacturing and other fields



**Horizontal Airflow
HCB-1600H**

- HEPA, with high efficiency up to 99.99% @ 0.3μm
Special design for conducting filter integrity testing
- Ergonomic design with top air inflow
Built-in interior lighting, to avoid eyestrain
- Special design for ensuring uniform air outflow
- Pre-cleaning function
- One-piece drug hanger made of stainless steel
- One-piece workbench made of 304 stainless steel 6mm ultraviolet-proof tempered glass
Special design for ensuring uniform airflow along back plate
- UV lamp interlocking
Patented technology UV lamp one-touch appointment/startup delay patented technology
- Universal caster for moving conveniently
Special underframe design for ensuring high stability

Product Advantages

-  Moisture-proof and flame-retardant glass fiber high efficiency particulate air filters (HEPA), ensure a filtration efficiency of 99.99% or above for particles of 0.3μm or above
-  With the industry's first pre-cleaning function: namely pre-cleaning the workspace before sample treatment so as to further improve the protection of samples
-  The air cleanliness conforms to Class 5 of the ISO14644.1 Standard, and is better than Class 100 cleanliness requirements
-  **One-key operation, convenient and secure:**
 - With a timed UV lamp turn off function, it is adjustable between 0 and 99 minutes
 - With the lighting lamp and UV lamp interlock function; the UV lamp can only be lit when the lighting lamp is off. The UV lamp is turned off immediately if the fluorescent lighting is turned on to prevent mis-operation
-  Built-in lighting prevents exposure to fluorescent lamp to avoid eye stain



**Vertical Airflow
HCB-900V/1300V**

- HEPA, with high efficiency up to 99.995% @ 0.3μm
Special design for conducting filter integrity testing
- Ergonomic design with top air inflow
Built-in interior lighting, to avoid eyestrain
- Special design for ensuring uniform air outflow
- Pre-cleaning function
- Ergonomic design of drop-down front sash
- One-piece workbench made of 304 stainless steel
6mm ultraviolet-proof tempered glass
- UV lamp interlocking
Patented technology UV lamp one-touch appointment/startup delay, patented technology
- Universal caster for moving conveniently
Heavy-duty stand

Laminar Flow Cabinet

Specifications

| Model | | HCB-900V | HCB-1300V | HCB-1600H | | |
|-----------------------|-----------------------------------|--|--|--|--|------------------|
| General Specification | Flow Type | Vertical | Vertical | Horizontal | | |
| | Power Supply(V/Hz) | 220/50 | 115/60 | 220/50/60 | | |
| | Power (W) | 1200 | 1200 | 350 | | |
| | Vibration Amplitude (UM) | 2 | 2 | 2 | | |
| | Exhaust Filter Typical Efficiency | H13 HEPA,99.99%@0.3um | H13 HEPA,99.99%@0.3um | H13 HEPA,99.99%@0.3um | | |
| | Average Velocity (M/S) | 0.2-0.4 | 0.2-0.4 | 0.2-0.4 | | |
| | Lighting Intensity (Lux) | ≥300 | ≥300 | ≥1000 | | |
| | Sound Level (dB(A)) | 60 | 60 | 61 | | |
| | Sash Opening(mm) | Max 310 | Max 310 | Max 310 | | |
| Dimensions | Net/Gross Weight | kg | 115/145 | 145/171 | 145/171 | 165/214 |
| | | lbs | 254/319 | 320/376 | 320/376 | 363.7/471 |
| | Internal Dimensions (W*D*H) | mm | 900*530*520 | 1300*530*520 | 1300*530*520 | 1710*550*750 |
| | | in | 35.4*20.9*20.5 | 51.2*20.9*20.5 | 51.2*20.9*20.5 | 67.3*21.7*29.6 |
| | External Dimension (W*D*H) | mm | 970*630*1730 | 1370*630*1730 | 1370*630*1730 | 1780*790*1960 |
| | | in | 38.2*24.8*68.1 | 53.9*24.8*68.1 | 53.9*24.8*68.1 | 70.1*31.1*77.2 |
| | Packing Dimensions (W*D*H) | mm | 1100*750*1290 | 1510*770*1280 | 1510*770*1280 | 1875*940*1410 |
| | | in | 43.3*29.5*50.8 | 59.4*30.3*50.4 | 59.4*30.3*50.4 | 73.8*37.0*55.5 |
| | Support Stand(mm) | | 755mm (standard) | 755mm (standard) | 755mm (standard) | 765mm (standard) |
| | Container Load (20'40'40'H) | | 15/33/33 | 10/22/22 | 10/22/22 | 6/12/12 |
| Accessories | UV Lamp | Y | Y | Y | Y | |
| | Power Outlet Socket | 2 Sockets | 2 Sockets | 2 Sockets | / | |
| | Conversion Socket | Optional | Optional | Optional | / | |
| | Cleanliness Classification | ISO 14644.1 Class 5 | ISO 14644.1 Class 5 | ISO 14644.1 Class 5 | ISO 14644.1 Class 5 | |
| | Valve Port | Optional | Optional | Optional | Optional | |
| | Side Wall Service Taps | Optional gas tap, vaccum tap, compressed air tap and water tap | Optional gas tap, vaccum tap, compressed air tap and water tap | Optional gas tap, vaccum tap, compressed air tap and water tap | Optional gas tap, vaccum tap, compressed air tap and water tap | |
| | Antimicrobial Coating | Optional | Optional | Optional | Optional | |
| | 316SS Work Surface | Optional | Optional | Optional | Optional | |
| | Work Surface | Single-piece | Single-piece | Single-piece | Single-piece | |
| | Touch Screen | / | / | / | / | |
| | Motorized Window | / | / | / | / | |
| | Foot Switch | / | / | / | / | |
| | Armrest | Y | Y | Y | Y | |
| | Exhaust Duct Connection Kits | / | / | / | / | |
| Others | Certification | CE | / | CE | CE | |



HCB-1300VS/1600VS

Product Advantages

- Positive pressure chamber design ensures air velocity uniformity and forms effective seal
- Ergonomic cabinet design with a 10° inclination, to improve user comfort during operation
- The workspace is made of easy-to-clean, high-quality, anti-corrosion 304 stainless steel
- The front sash window and the left & right-side glass panels are all non-reflective tempered UV glass to maximise natural light during operation, which also makes the glass easy to clean
- Designed with access ports on left and right sides, easy and convenient to mount valves such as gas supply

*Haier Biomedical reserves the right to change products and specifications without prior notice

Laminar Flow Cabinet

Specifications

| Model | | HCB-1300VS | HCB-1600VS | |
|-----------------------------|-----------------------------------|--|--|----------------|
| General Specification | Flow Type | Vertical | Vertical | |
| | Power Supply(V/Hz) | 220/50/60 | 220/50/60 | |
| | Power (W) | 1350 | 1200/1350 | |
| | Vibration Amplitude (UM) | 2 | 2 | |
| | Exhaust Filter Typical Efficiency | H13 HEPA,99.99%@0.3um | H13 HEPA,99.99%@0.3um | |
| | Average Velocity (M/S) | 0.20-0.50 | 0.20-0.50 | |
| | Lighting Intensity (Lux) | ≥600 | ≥1000 | |
| | Sound Level (dB(A)) | 60 | 63 | |
| Dimensions | Sash Opening(mm) | Max 360 | Max 360 | |
| | Net/Gross Weight | kg | 180/232 | 202/260 |
| | | lbs | 397/511.5 | 445.3/573 |
| | Internal Dimensions (W*D*H) | mm | 1300*510*550 | 1600*510*550 |
| | | in | 51.2*20.1*21.7 | 60.0*20.1*21.7 |
| | External Dimesion (W*D*H) | mm | 1370*640*1820 | 1670*640*1820 |
| | | in | 53.9*25.2*71.7 | 65.7*25.2*71.7 |
| | Packing Dimensions (W*D*H) | mm | 1480*845*1290 | 1780*840*1270 |
| | | in | 58.3*33.3*50.8 | 70.1*33.1*50.0 |
| | Support Stand(mm) | 750mm (standard) | 750mm (standard) | |
| Container Load (20'40'40'H) | 11/20/20 | 7/14/14 | | |
| Accessories | UV Lamp | Y | Y | |
| | Power Outlet Socket | 2 sockets | 2 sockets | |
| | Conversion Socket | Optional | Optional | |
| | Cleanliness Classification | ISO 14644.1 Class 4 | ISO 14644.1 Class 4 | |
| | Valve Port | 2 on each side | 2 on each side | |
| | Side Wall Service Taps | Optional gas tap, vaccum tap, compressed air tap and water tap | Optional gas tap, vaccum tap, compressed air tap and water tap | |
| | Antimicrobial Coating | Optional | Optional | |
| | 316SS Work Surface | Optional | Optional | |
| | Work Surface | Single-piece | Single-piece | |
| | Touch Screen | / | / | |
| Others | Motorized Window | / | / | |
| | Foot Switch | / | / | |
| | Armrest | Y | Y | |
| | Exhaust Duct Connection Kits | / | / | |
| | Certification | / | / | |

*Haier Biomedical reserves the right to change products and specifications without prior notice

Ductless Fume Hood

Scope of Application

Haier Biomedical Ductless Fume Hood utilizes carbon filtration to protect laboratory personnel and the environment from toxic chemical fumes, odours, and particles. It is widely used in universities, public health, chemical and petrochemical industry, drug testing, museums and many other fields which involve the use of chemicals

10-inch touch screen

Real-time display of VOC concentration, wind speed, filter service life and UV lamp running time. Alerts when filter or UV lamp lifespan falls below 10%

EC fan

Energy-efficient EC fan for a longer service life and reliable and quieter operation

Air velocity sensor

Intelligent constant air velocity ensures the stability


PIR detection

In intelligent mode, the PIR sensor will detect when the operating area has been clear of personnel for more than 15 minutes and automatically switch to LNS green-saving mode

Universal casters with built-in threaded support legs

Electrical sockets

Included as standard are two, waterproof sockets with timing technology, which allows the user to program timed on/off function



HTF-1200W

Optional Components



UV lamp



Dual surveillance cameras




Electrically operated glass sash



IoT module




Height adjustable stand





PVC worktable


Ductless Fume Hood


Product Advantages

 Intelligent constant speed design, ensures the best adsorption effect of the filter


 Standard VOC sensor detects presence of volatile organic compounds in the exhaust and alarms to indicate filter saturation

 Energy-efficient EC fan for a longer service life and reliable operation

 Multiple stackable filter types: remove vapor contaminants from the air, can be combined to handle different types of chemicals and particulates during the same application

 Offers a broad range of accessories and options to meet clients' different requirements

Product Dimension



Types of Filters

Carbon filters are used to remove toxic gases, hazardous fumes, and odors. These filters are constructed from high-quality carbon pellets and durable chemical-resistant cases. Below is Haier Biomedical filters list

| Code No | Name | Suitable Applications |
|---------|---|---|
| A | Standard Filter (Mixed with N4A1, N4B1 and N4G1 activated carbon) | Suitable for all common laboratory chemicals, especially with organics. When no specific requirements are present, or when more than one type of chemical is used |
| B | Organic Solvent Filter (N4G1 activated carbon) | Removal of toluene, benzene, xylene, acetone, acetic acid, carbon tetrachloride, chloroform, CxHy, VOC, etc |
| C | Acid Gas Filter (N4A1 activated carbon) | Removal of hydrochloric acid, sulfuric acid, hydrofluoric acid, hydrogen sulfide, sulfur dioxide, etc |
| D | Alkaline Gas Filter (N4B1 activated carbon) | Removal of ammonia, amine, etc. |
| E | Formaldehyde Filter | Mainly used for adsorption of formaldehyde |
| F | Particle Filter | Mainly used to absorb dust |

Specifications

| Model | HTF-1200W | |
|-----------------------------------|--------------------------|--|
| Power Supply (V/Hz) | 220/50/60 | |
| Power (W) | 1200 | |
| Sound Level (dB(A)) | ≤58 | |
| Inflow Air Velocity (m/s) | 0.4-0.6 | |
| Fluorescent Light Intensity (lux) | ≥1000 | |
| Filtration | First Filter | Activated Carbon with Granular Media bed (6 different filter types available, codes A-F) |
| Elements | Second Filter (optional) | Activated Carbon with Granular Media bed (6 different filter types available, codes A-F) |
| Exterior Dimensions (W*D*H)(mm) | 1336*780*2120 | |
| Operation Dimensions (W*D*H)(mm) | 1230*760*655 | |
| Packaging Dimensions (W*D*H)(mm) | 1400*925*1665 | |
| Net Weight (Kg) | 230 | |
| Packaging Weight (Kg) | 290 | |
| Display | 10 inch screen | |
| VOC Sensor | Standard (1) | |
| Certification | CE | |

*Haier Biomedical reserves the right to change products and specifications without prior notice.