



Medical Oxygen Concentrator

Scope of Application:

Medical oxygen administered by the concentrator is beneficial to cure respiratory, heart, cardiovascular and chronic pulmonary diseases, as well as other diseases requiring concentrated oxygen as key examples.

Innovative Design

- Low maintenance
- Easy to maintain
- Portable and user friendly

Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China Tel: +86-0532-88935593 Website: www.haiermedical.com













Product Design and Advantages 〈







Dual Flow Model

High quality humidifier with 6PSI pressure release



Multiple Alarms

Multiple alarms are optional to guarantee stable performance



Innovative "Cooling"

Safe and secure, improved reliability, extended life of unit



Sieve Bed

Imported lithium molecular sieve durable, high efficiency



Large LCD Display

Operates continuously for 24h* 365days



Medical Grade Material

Environmentally friendly material



Specifications



Model	JAY-10W
Flow Rate (LPM)	0-10
Purity	93% (±3%)
Outlet Pressure (MPA)	0.04-0.07 (6-10PSI)
Sound Level Power (dB (A))	≤50
Consumption (W)	≤880
LCD Display	On/Off, operating pressure, current working time, accumulated working time, preset timer from 10mins to 40hours
Alarm	Power failure alarm
Net Weight (kg)	27
Size (mm)	365×375×600

(Options: temperature alarm, low purity alarm, extra humidifier, pulse oximeter and nebulizer.)