



# Medical Oxygen Concentrator

## Scope of Application:

Medical oxygen administered by the concentrator is beneficial to cure respiratory, heart, cardiovascular and chronic pulmonary diseases, as well as other diseases requiring concentrated oxygen as key examples.

# Innovative Design

- Low maintenance
- Easy to maintain
- Portable and user friendly

#### Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China Tel: +86-0532-88935593 Website: www.haiermedical.com













## Product Design and Advantages







Higher outlet pressure to ensure optimum oxygen flow delivery



Low maintenance, easy to maintain, portable and user friendly



Optional pulse oximeter for measuring the oxygen saturation

# Specifications (

Model	JAY-5BW
Flow Rate (LPM)	0-5
Purity	93% (±3%)
Outlet Pressure (MPA)	0.04-0.07 (6-10PSI)
Sound Level Power (dB (A))	≤38
Consumption (W)	≤350
LCD Display	On/Off, operating pressure, current working time, accumulated working time, preset timer from 10mins to 40hours
Alarm	Power failure alarm
Net Weight (kg)	16
Size (mm)	360×243×600