

# Smart Mobile Clinic Inoculation Solution



## Convenient inoculation clinic, mobile health care security

The health care solution improving medical services by reaching remote and hard to reach locations, such as towns, small villages and mountainous regions.

**Qingdao Haier Biomedical Co.,Ltd.**

No.280 Feng Yuan Road, High-tech Zone,  
Qingdao, 266109, P.R. China  
Tel: +86-0532-88935593  
Website: [www.haiermedical.com](http://www.haiermedical.com)



Haier Biomedical  
International



Haier Biomedical  
International



@haiermedicalint



@haierbiomedicalint



Haier Biomedical  
International

# Smart Mobile Clinic Inoculation Solution

## Social Benefits



### Improved and convenient vaccination services

Innovative and modern solution service for mobile vaccination needs, the Haier Biomedical comprehensive mobile solution is fitted with critical medical equipment generally found in a clinic.



### Promote and improve the equity in immunization services

Improve the quality of health services for people in remote, rural and disaster-hit areas, and ensure that children in these regions as well as migrant children receive the same medical treatment as urban children, thus promoting fairness in immunization services.



### Improve the access to immunization services


The mobile vaccination platform has extended its coverage to more remote mountainous areas and villages.



### Guarantee the accuracy and immediate reporting of data and information remotely

Digital inoculation management. Records of the whole process are traceable. The inoculation data is connected with the disease control system in real time to ensure the accuracy and timeliness of the information.

## Solution



## Parameters

Item	Inoculation Refrigerator			Solar Direct Drive Combined Refrigerator/Freezer			
Performance	Operating Temp (°C)	Storage Temp (°C)	Vaccine Capacity (vial)	Vaccine Capacity (vial)	Refrigerating Volume (L)	Freeze Protection Grade	Freezing Volume (L)
Parameter	5~43	2~8	200	300	37.5	A	32

Item	Power Supply (battery)		Medical Vehicle						
Performance	Duration	Charging Time (h)	Double Tanks (L)	Fuel Consumption (L/km)	Maximum Mileage (km)	Length (mm)	Width (mm)	Height (mm)	Weight (t)
Parameter	3 refrigerators running for 1.5days (6h/day)	10	2*100	18/100	100	8995	2500	3680	11.9