

The New Force in Centrifugation

Haier Biomedical research centrifuge packages help laboratories process more samples with confidence. Each package combines the centrifuge, rotor, buckets and adapters required for common research workflows, from cell culture and conical tube processing to high-speed microcentrifugation.

With intuitive controls, high-speed performance and refrigeration options, the range supports repeatable results, efficient workflows and sample protection in demanding modern research laboratories.

Secure 2026 package pricing today.

Call 01932 780 070 or email
sales@haierbiomedical.co.uk



Floor Standing Centrifuge



Part Code	LX-210L750R Floor-Standing Refrigerated Centrifuge Package	Package Price*
CentpackFloor1    	LX-210L750 Floor Standing Centrifuge Pack. Includes: LX-210L750R Floor Standing Refrigerated Centrifuge Fixed-angle rotor 8 x 50ml Fixed-angle rotor 6 x 250ml Fixed-angle rotor 6 x 500ml	£39,484.00 £19,999.00 Almost Half Price!

The Haier Biomedical LX-210L750R floor-standing refrigerated centrifuge package is built for demanding research laboratories that need high-speed, large-volume fixed-angle processing. Supplied with 8 x 50 ml, 6 x 250 ml and 6 x 500 ml rotors, it gives teams the flexibility to support multiple workflows from one platform.

- Tool-less Rotor exchange for safe and fast switching of rotors
- Includes three fixed-angle rotors for 50 ml, 250 ml and 500 ml processing
- Powerful and fast refrigeration helps protect temperature-sensitive samples during runs
- 7" touchscreen with named programs and RFID rotor recognition for consistent operation

*Price valid until 31st December 2026, excludes VAT and delivery.



Cell Culture Packages



Part Code	LX-155T750R Refrigerated Cell Culture Package	Package Price*
CentpackRCC3	LX-155T750R 4 x 750 Refrigerated Cell Culture Pack. Includes: LX-155T750R Refrigerated Centrifuge QX-750 Swing-out Rotor QX-750 500ml Buckets (4) Adapters for 28 x 50ml Conical Tubes Adapters for 64 x 15ml Conical Tubes	£11,827.00 £7,096.20 Save 40%

The Haier Biomedical LX-155T750R refrigerated cell culture package gives research laboratories the capacity and temperature control needed for high-volume cell culture, media preparation and routine sample processing. With a 4 x 750 ml swing-out rotor and conical tube adapters, it flexibly processes 50 ml and 15 ml tubes while helping protect sample integrity in temperature-sensitive workflows.

- Processes up to 28 x 50 ml or 64 x 15 ml conical tubes
- 7" touchscreen with named programmes and RFID rotor recognition for safe, repeatable operation
- ECO fastBlast refrigeration pre-cools in as little as 5 minutes, reducing wait time and energy use
- High-speed processing up to 15,500 rpm with optional fixed-angle rotors



Part Code	LX-60T750R Refrigerated Cell Culture Package	Package Price*
CentpackRCC4	LX-60T750R 4 X 750 Refrigerated Cell Culture Pack. Includes: LX-60T750R Refrigerated Centrifuge QX-750 Swing-out Rotor QX-750 500ml Buckets (4) Adapters for 28 x 50ml Conical Tubes Adapters for 64x 15ml Conical Tubes	£10,907.00 £6,544.20 Save over 40%

The Haier Biomedical LX-60T750R refrigerated cell culture package is a cost-effective, high-capacity solution for routine research laboratories. With a 4 x 750 ml swing-out rotor, buckets and adapters for common conical tubes, it delivers reliable separation, simple operation and excellent tube capacity in a compact benchtop format.

- Swing-out rotor package for routine cell culture and conical tube processing
- Simple, high-contrast touchscreen for fast parameter changes
- Compact, low-height design supports easier loading and efficient use of bench space



Part Code	LX-155T500R Refrigerated Cell Culture Package	Package Price*
CentpackRCC2	LX-155T500R 4 x 500 Refrigerated Cell Culture Pack. Includes: LX-155T500R Refrigerated Centrifuge QX-500 Swing-out Rotor QX-500 500ml Buckets (4) Adapters for 16 x 50ml Conical Tubes Adapters for 40 x 15ml Conical Tubes	£10,575.00 £6,345.00 Save 40%

The Haier Biomedical LX-155T500R refrigerated centrifuge is designed for research laboratories that require powerful temperature control, high-capacity processing and advanced safety features. It helps protect sample integrity, improve workflow efficiency and support reliable scientific results.

- 7" touchscreen with named programmes and RFID rotor recognition
- ECO fastBlast refrigeration pre-cools in as little as 5 minutes, saving time and energy
- 14 additional rotors available separately for greater application flexibility
- Low noise level of <55 dB(A) to support a quieter working environment



Part Code	LX-60T500 Ventilated Cell Culture Package	Package Price*
CentpackVCC2	LX-60T500 4 X 500 Ventilated Cell Culture Pack. Includes: LX-60T500 Ventilated Centrifuge QX-500 Swing-out Rotor QX-500 500ml Buckets (4) Adapters for 16 x 50ml Conical Tubes Adapters for 40 x 15ml Conical Tubes	£6,775.00 £4,065.00 Save 40%

The Haier Biomedical LX-60T500 ventilated centrifuge is designed for routine research laboratories that need reliable separation, excellent tube capacity and strong value. Its simple interface and compact footprint make it ideal for everyday cell culture and conical tube workflows.

- Simple, intuitive, high-contrast LCD display with buttons and control knob
- Low-height, compact design helps save valuable bench space
- Ventilated design supports efficient heat dissipation during centrifugation



Part Code	Description	Package Price*
CentpackRCC1	LX-185T100R 4 x 100 Refrigerated Cell Culture Pack. Includes: LX-185T500R Desktop Refrigerated High-speed Centrifuge QX-100 Swing-out Rotor QX-100B 100ml Buckets (4) Adapters for 4 x 50ml Conical Tubes Adapters for 8 x 15ml Conical Tubes	£7,791.00 £4,674.60 Save 40%

The Haier Biomedical LX185T100R, a state-of-the-art refrigerated centrifuge for research labs, offers precise temperature control, high-capacity processing, and advanced safety features to ensure sample integrity, optimize efficiency, and support scientific discoveries.

- Powerful refrigeration, ECOfastblastpre-cool in as little as 5 minutes, saving time and energy
- Low Noise Level of <52dB(A) to ensure a great working environment
- High Speed Performance – up to 18,500 rpm / 30,130xg (with optional rotors)
- 7" Touchscreen with named programs and RFID Rotor tracking to ensure safety



Part Code	Description	Package Price*
CentpackVCC1	LX-60T100 4 X 100 Ventilated Cell Culture Pack. Includes: LX-60T100-J Ventilated Centrifuge QX-100 Swing-out Rotor QX-100M Buckets 4 x 50ml Conical Tubes QX-100M Buckets 8 x 15ml Conical Tubes	£3,546.00 £2,127.60 Save 40%

The Haier Biomedical LX60T100 is a ventilated centrifuge for research and diagnostic labs, offering efficient heat dissipation, precise speed control, and a compact design. It's an ideal choice for labs needing high performance in a space-saving solution.

- High-contrast screen with easy-to-use buttons and a knob for quick and effortless parameter adjustments.
- Compact design to optimise lab space and a low-profile design.
- Low Noise Level of <56dB(A) to ensure a great working environment

Microcentrifuge Packages



Part Code	LX-165T2R Refrigerated High-Speed Microcentrifuge Package	Package Price*
CentpackRM1	Refrigerated High-Speed Microcentrifuge Pack. Includes: LX-165T2R High-Speed Refrigerated Micro centrifuge QFA-26 24 x 2ml Fixed Angle Rotor with Click & Spin Lid	£5,288.00 £3,172.80 Save 40%
CentpackRM1B	Refrigerated High-Speed Microcentrifuge Pack. Includes: LX-165T2R High-Speed Refrigerated Micro centrifuge QFA-26 24 x 2ml Fixed Angle Biocontainment Rotor	£5,438.00 £3,262.80 Save 40%

The Haier Biomedical LX165T2R is the flagship refrigerated high-speed microcentrifuge, ideal for molecular biology research. Advanced temperature control and user-friendly operation ensure the highest level of precision and reliability for sensitive sample handling.

- High-speed performance up to 16,500 rpm / 25,872 x g, rotor dependent
- 7" touchscreen with named programmes and RFID rotor tracking for safer operation
- ECO fastBlast refrigeration pre-cools in just 5 minutes, saving time and energy
- Soft-close lid for smooth, quiet closing



Part Code	LX-165T2-J Ventilated High-Speed Microcentrifuge Package	Package Price*
CentpackVM1	Ventilated High-Speed Microcentrifuge Pack. Includes: LX-165T2-J High-Speed Ventilated Microcentrifuge QFA-26 24 x 2ml Fixed Angle Rotor with Click & Spin Lid	£2,088.00 £1,252.80 Save 40%
CentpackVM1B	Ventilated High-Speed Microcentrifuge Pack. Includes: LX-165T2-J High-Speed Ventilated Microcentrifuge QFA-26 24 x 2ml Fixed Angle Biocontainment Rotor	£2,038.00 £1,222.80 Save over 40%

The Haier Biomedical LX-165T2-J ventilated high-speed microcentrifuge is designed for molecular biology and low-volume sample processing. Efficient ventilation, intuitive controls and robust construction support consistent performance and precise results.

- High-speed performance up to 16,500 rpm / 25,872 x g, rotor dependent
- LED keypad and control knob for simple, intuitive operation
- Wide viewing-angle display for easier monitoring
- High-quality construction, including a premium 316 stainless steel bowl

Did you know that Haier Biomedical can service a wide range of laboratory equipment, including centrifuges?

At Haier Biomedical, we understand that your laboratory equipment plays a crucial role in the success of your work. That's why we offer a comprehensive range of maintenance services designed to keep your equipment in peak condition, ensuring maximum uptime, compliance, and long-term reliability.

Service Plan Tailored to Your Needs - Choose Your Plan

Recommended option for units still within warranty

Part Code	GOLD Fully Comprehensive	SILVER Unlimited Breakdown	BRONZE Preventative Maintenance
One preventative maintenance visit PM Checklist includes: Preventative maintenance, critical component review and operating environment checks	Yes	Yes	Yes
Free consumables Replacement batteries and filters*	Yes	Yes	Yes
Travel	Inclusive	Inclusive	-15%
Labour	Inclusive	Inclusive	-15%
Parts	Inclusive	-10%	-5%
Priority call outs	Yes	Yes	No

*Applicable to ULT freezers only, subject to make and model. Replacement during service visit.

Haier Biomedical UK offers three levels of service contracts tailored to suit your individual service and maintenance needs. We have a dedicated UK based service support team and factory trained engineers.

Haier Biomedical UK can also provide flexible bespoke contracts and include services additional preventative maintenance visits, annual calibrations, temperature mapping and IQ/OQ/PQ services.

Products Covered

- Cryogenic Freezers (-150°C)
 - ULT Freezers (-80°C)
 - Freezers (-20°C to -40°C)
 - Refrigerators (+2°C to +8°C)
 - Blood Refrigerators
 - Plasma Freezers
- Liquid Nitrogen Storage (-196°C)
 - Cold and Freezer Rooms
 - Centrifuges
 - CO₂ Incubators
 - Standard Incubators
 - Climate Chambers
- Class II Safety Cabinets
 - Laminar Flow Cabinets
 - Fume Hoods
 - Autoclaves
 - Sterilisers
 - Temperature Monitoring

